



MANUFACTURER AND DISTRIBUTOR OF MEDICAL EQUIPMENT

PRODUCT CATALOG



*We don't see limitations in a man – we see efficiency
restored through modern solutions*



■ QUALITY

With full commitment, we are constantly setting new and better technical standards, which have made a significant contribution to improving the lives of those who use the company's products.

Reha Fund products meet both Polish and EU legal requirements for medical equipment. Furthermore, we can boast about the implementation of the ISO 13485 standard, which is an international quality standard for medical devices.

The certification for compliance with this standard confirms that our company has a management system that ensures the high quality of our products, the necessary resources for maintaining and improving the system, as well as complies with all legal requirements and is ready to implement the new MDR regulations, which explicitly oblige us to maintain an effective quality management system.

Through the Quality Management System, we ensure that the requirements of our customers are 100% satisfied.





We have been operating on the Polish and European medical market since 1992. We are a manufacturer, distributor and importer of medical devices and, as one of the first Polish companies, we have shaped the Polish market by introducing both our own innovative products and those of world leaders in the industry.

We believe that the crucial thing is the satisfaction of our customers. Therefore, we place great importance not only on the quality of our products and services, but also on the daily service of each customer. Such an approach ensured us the trust of our many satisfied customers.

■ MISSION

Together with our business partners, we build a system of services, patient care and distribution of specialized medical products.

We do not see any limitations in a man – we see efficiency restored through modern solutions.

■ VISION

We consistently strive to be one of the market leaders in rehabilitation and medical products increasing the quality and comfort of life.

■ VALUES

Satisfaction – better and more efficient life of our customers.

Responsibility – thanks to the investment in new solutions and products, we feel responsible for the health of thousands of patients.

Respect – the ability to listen and understand our customers, business partners as well as company's employees.

Flexibility – an individual approach to the needs of our customers, employees and business partners.

Trust – we have been present on the market of rehabilitation and medical products for decades. Thousands of patients and institutions have already trusted us. We believe in long-term ventures.

TABLE OF CONTENTS

OXYGEN THERAPY	7
PORTABLE OXYGEN CONCENTRATORS	9
STATIONARY OXYGEN CONCENTRATORS	10
ACCESSORIES FOR STATIONARY OXYGEN CONCENTRATORS	12
FINGERTIP PULSE OXIMETERS	14
MEDICAL SUCTION UNITS	15
SLEEP APNEA	17
CPAP MACHINES	19
MASKS	23
ORTHOPEDIC PRODUCTS	27
FOOT AND ANKLE	29
KNEE	33
ABDOMEN	38
SPINE	40
SHOULDER	44
ELBOW	46
WRIST AND HAND	48
NECK	53
REHABILITATION PRODUCTS	55
AUXILIARY AND SANITARY AIDS	63
SANITARY AIDS	65
AUXILIARY AIDS	79
PILL ORGANIZERS AND ACCESSORIES	81

ANTI-BEDSORE MATTRESSES	87
ANTI-BEDSORE PNEUMATIC MATTRESSES.....	89
ANTI-BEDSORE FOAM MATTRESSES	95
COVERS	98
UNDERPADS	99
ORTHOPEDIC PILLOWS FOR SLEEPING	101
ANTI-BEDSORE SEAT CUSHIONS	107
ANTI-BEDSORE PNEUMATIC SEAT CUSHIONS	109
ANTI-BEDSORE FOAM SEAT CUSHIONS	111
ANTI-BEDSORE AND POSITIONING PADS	113
WALKING AIDS.....	121
WALKING FRAMES, FOREARM WALKING FRAMES, HANDCARTS.....	123
4-WHEELED WALKERS	125
3-WHEELED WALKERS	128
CRUTCHES.....	129
TRIPOD, TETRAPOD, CANES	130
REHABILITATION EXERCISERS.....	131
WHEELCHAIRS.....	133
MOBILE LIFTS.....	146





OXYGEN THERAPY



Long-trusted distributor and authorized service of oxygen concentrators

We are a long-term partner of Drive DeVilbiss, Aerti and Kingon companies, and their products have been present in our offer for several years. Our vast experience not only in selling, but also in servicing oxygen concentrators has shown us one thing – these devices are one of the best in the world. Their very simple operation, high quality and good prices are the answer to all the patient's needs.

Rental of oxygen concentrators

Oxygen concentrators are essential for people with low blood oxygen levels. The purchase of such a device is associated with a large, one-time expense, as well as with the need to provide technical and service care. Therefore, renting an oxygen concentrator offered by the Reha Fund rental company is an ideal and economical solution for people who need quick access to oxygen.

We are a manufacturer and distributor of medical equipment, and for many years we have been renting oxygen concentrators and other medical equipment in Warsaw and all over Poland. Out of concern for the comfort of life and health of our customers, we **rent** only **modern oxygen concentrators** that are very easy to use, quiet and reliable. Due to their solidity we can offer a low price for their rental compared to the quality.

Contact us

Our specialists will help you choose the best device.
Shipment is carried out on the next business day!

Find out more!

Visit our website www.rehafund.pl

PORTABLE OXYGEN CONCENTRATOR



iGO2

Available colors:

● gray

- Lightweight
- With a SmartDose function that automatically adjusts the amount of administered oxygen depending on current needs
- With a triple power supply system: AC, DC and an internal rechargeable battery
- Produces highly concentrated oxygen over the entire flow range

TECHNICAL PARAMETERS	
dimensions	213 × 89 × 218 mm
noise level	< 37 dB
weight	2.2 kg
oxygen concentration	90% (+4/-3%) (at each flow)
adjustable oxygen flow	1 to 5 l/min for standard flow, 1 to 4 l/min for SmartDose mode flow
DC power supply	13.8 – 17 VDC (internal battery, car cigarette lighter)
AC power supply	90-240 VAC, 50-60 Hz
internal battery run time	3.5 h (at 2 l/min flow)

Intended use:

Intended for people with reduced oxygen content in the blood, especially all patients with COPD, lung cancer, COVID-19. Intended for people who need constant oxygen support from the outside, but while being able to move freely. This device makes it possible to move outside the home.



PORTABLE OXYGEN CONCENTRATOR

kingon

KINGON P2

Available colors:

○ white

- Lightweight
- With a triple power supply system: AC, DC and an internal rechargeable battery
- Produces highly concentrated oxygen over the entire flow range
- Works only in a pulse mode

TECHNICAL PARAMETERS	
dimensions	221 × 85 × 160 mm
noise level	48 dB
weight	1.98 kg
oxygen concentration	90% (-3%/+6%)
adjustable oxygen flow	1 to 5 l/min
DC power supply	14.4 VDC (internal battery, car cigarette lighter)
AC power supply	90-240 VAC, 50-60 Hz
internal battery run time	up to 4 h

Intended use:

The Kingon P2 oxygen concentrator is used as prescribed by a physician for patients requiring supplemental oxygen. It delivers highly concentrated oxygen and is used with a nasal cannula to transport the oxygen from the concentrator to the patient. Kingon P2 is small, portable and can be used at home, in health centers and in various mobile environments.



OXYGEN CONCENTRATOR – DeVilbiss



525KS

- Quiet
- Lightweight
- Maneuverable
- Produces highly concentrated oxygen

TECHNICAL PARAMETERS	
dimensions	62.2 × 34 × 30.4 cm
noise level	< 40 dB
weight	16.3 kg
oxygen concentration	93% +/-3%
adjustable oxygen flow	0.5 – 5.0 l/min
power consumption	290 W at 2 l/min flow, 312 W at 5 l/min flow

Intended use:

Medical use: hospital and home use, medical units in areas where oxygen is difficult to deliver, assisting in the treatment of lung disease, obstructive pulmonary disease, pneumothorax and air starvation. Contraindications: disabled patients with severe carbon monoxide poisoning.



Available colors:

● black

OXYGEN CONCENTRATOR – DeVilbiss



1025KS

- High oxygen flow
- Maneuverable
- Produces highly concentrated oxygen

TECHNICAL PARAMETERS	
dimensions	62.2 × 34 × 30.4 cm
noise level	< 57 dB
weight	19 kg
oxygen concentration	93% +/-3%/-6%
adjustable oxygen flow	2.0 – 10.0 l/min
power consumption	664 W – average consumption

Intended use:

Medical use: hospital and home use, medical units in areas where oxygen is difficult to deliver, assisting in the treatment of lung disease, obstructive pulmonary disease, pneumothorax and air starvation. For patients requiring increased oxygen ventilation. Contraindications: disabled patients with severe carbon monoxide poisoning.



Available colors:

● black

OXYGEN CONCENTRATOR – AERTI



AE-5W

Available colors:

● dark gray

- Quiet
- Maneuverable
- Produces highly concentrated oxygen

TECHNICAL PARAMETERS	
dimensions	314 × 240 × 515 mm
noise level	< 36 dB
weight	21.5 kg
oxygen concentration	93% +/-3%
adjustable oxygen flow	0.5 – 5.0 l/min
power consumption	350 W

Intended use:

Medical use: hospital and home use, medical units in areas where oxygen is difficult to deliver, assisting in the treatment of lung disease, obstructive pulmonary disease, pneumothorax and air starvation. For patients requiring increased oxygen ventilation. Contraindications: disabled patients with severe carbon monoxide poisoning.



OXYGEN CONCENTRATOR – AERTI



AR-5

Available colors:

○ white

- Quiet
- Maneuverable
- Produces highly concentrated oxygen

TECHNICAL PARAMETERS	
dimensions	314 × 240 × 515 mm
noise level	< 46 dB
weight	13.9 kg
oxygen concentration	93% +/-3%
adjustable oxygen flow	0.5 – 5.0 l/min
power consumption	350 W

Intended use:

Medical use: hospital and home use, medical units in areas where oxygen is difficult to deliver, to support the treatment of lung disease, obstructive lung disease, pneumothorax and air starvation. For patients requiring increased oxygen ventilation. Contraindications: disabled patients with severe carbon monoxide poisoning.



Ask our Sales Representative about the additional equipment for stationary oxygen concentrators.

OXYGEN TUBING

02-1302-E – 2.1 m

2025 – 7.6 m

2050/1 – 10 m



OXYGEN MUSTACHE AND OXYGEN TUBING CONNECTOR

1215



OXYGEN MASK FOR ADULTS OR CHILDREN WITH OXYGEN TUBING 2.1 m

8115-0-50/1 – mask for adults

8115-0-50/2 – mask for children



OXYGEN MUSTACHE

1600 – mustache 2.1 m

1600/1 – mustache 5 m



HUMIDIFIER AND OXYGEN CONCENTRATOR ADAPTER

1215/1



SHORT TUBE FOR CONNECTING HUMIDIFIER TO OXYGEN CONCENTRATOR

E505.7



HUMIDIFIER FOR OXYGEN CONCENTRATOR

7600



FINGERTIP PULSE OXIMETERS

YK-81C

Available
colors:

● black

- Easy to use
- OLED display
- Non-invasive blood oxygen level monitoring
- Equipped with alarms to monitor abnormal saturation levels and heart rate

TECHNICAL PARAMETERS	
dimensions	58 × 36 × 33 mm
saturation measurement range	70-99% ($\pm 2\%$ at $>70\%$)
heart rate measurement range	30-240 BPM ($\pm 1\%$ at >100 BPM)
power supply	2 AAA batteries
battery run time	up to 40 h

Intended use:

This fingertip pulse oximeter is an innovative medical device with non-invasive and on-going SPO2 and pulse detection functions. Being a portable device, it is able to measure SPO2 and heart rate values quickly and accurately.



FINGERTIP PULSE OXIMETERS

K-1

Available
colors:

● yellow

- Easy to use
- OLED display
- Non-invasive blood oxygen level monitoring
- Equipped with alarms to monitor abnormal saturation levels and heart rate

TECHNICAL PARAMETERS	
dimensions	44 × 28 × 26 mm
saturation measurement range	70-99% ($\pm 2\%$ at $>70\%$)
heart rate measurement range	30-240 BPM ($\pm 1\%$ at >100 BPM)
power supply	internal rechargeable battery powered by USB connector
battery run time	working time on a charged battery up to 6 h

Intended use:

This fingertip pulse oximeter is an innovative medical device with non-invasive and on-going SPO2 and pulse detection functions. Being a portable device, it is able to measure SPO2 and heart rate values quickly and accurately.



MEDICAL SUCTION UNIT

7314P-NE



Available colors:

● blue

- Portable
- With a triple power supply system: AC, DC and an internal rechargeable battery
- High vacuum range adjustment

TECHNICAL PARAMETERS	
dimensions	21 × 20 × 21 cm
noise level	< 55 dB
maximum air flow	27 l/min
weight	3.0 kg
suction pressure range	50-550 mmHg (max. 780 kPa)
AC power supply	AC 100/240V, 50/60 Hz
internal battery run time	up to 4 h

Intended use:

Intended for people who require suction of body fluids, mucus and other discharge through the oral cavity, nose or larynx. It has a wide range of suction power, so it can be used by both adults and infants. It includes an internal (rechargeable) battery and a car cigarette lighter socket which makes transport and walking easier.



MEDICAL SUCTION UNIT

7E-C

yuwell

Available colors:

○ white

- Lightweight
- Adjustable suction power
- For domestic use

TECHNICAL PARAMETERS	
dimensions	370 × 180 × 280 mm
noise level	< 65 dB
maximum air flow	15 l/min
weight	4.4 kg
suction pressure range	200-760 mmHg
AC power supply	AC 100/240V, 50/60 Hz

Intended use:

A perfect household device for people who require suction of body fluids, mucus and secretions through the oral cavity, nose or larynx. The Yuwell company is a world-renowned recognized manufacturer of electromedical devices. High quality available at an affordable price.







SLEEP APNEA



Reha Fund is the official distributor of the most renowned devices for the treatment of sleep apnea

Our offer includes products from the world's best manufacturers, such as **Resmed** or **Yuwell**, helping in the conscious choice of the best device adjusted to the patient's needs.

We select and sell nasal and face masks, which are an indispensable equipment component in the treatment of sleep apnea.

Sleep apnea is extremely dangerous for the health and even life of the patient. It occurs in people who, as a result of excessive relaxation of the tissues of the upper respiratory tract, have limited or complete lack of airflow through the respiratory tract. People diagnosed with sleep apnea should get a CPAP machine, which, by providing a stream of air, clears the respiratory tract and eliminates the phenomenon of apnea.

Find out more!

Visit our website www.rehafund.pl

CPAP MACHINE – YUWELL

YH-580/YH-560

yuwell

Available colors:

- black
- white

- Automatic
- Detachable humidifier
- Quiet

TECHNICAL PARAMETERS	
dimensions	270 × 135 × 100 mm
noise level	30 dB
pressure range	4-20 mBar
weight	1.6 kg
RAMP function	Adjustable in the range 0-45 min
AC power supply	AC 100/240V, 50/60 Hz

Intended use:

The device is intended to deliver continuous positive airway pressure (CPAP) to the respiratory tract for the treatment of obstructive sleep apnea. An integrated, heated humidifier is used to increase air humidity in the CPAP machine, thus alleviating the symptoms of a dry nose or throat that may occur during the therapy.



CPAP MACHINE – YUWELL

YH-360

yuwell

Available colors:

- white

- Constant airway pressure
- Detachable humidifier
- Quiet

TECHNICAL PARAMETERS	
dimensions	270 × 135 × 100 mm
noise level	30 dB
pressure range	4-20 mBar
weight	1.6 kg
RAMP function	Adjustable in the range 0-45 min
AC power supply	AC 100/240V, 50/60 Hz

Intended use:

The device is intended to deliver continuous positive airway pressure (CPAP) to the respiratory tract for the treatment of obstructive sleep apnea. An integrated, heated humidifier is used to increase air humidity in the CPAP machine, thus alleviating the symptoms of a dry nose or throat that may occur during the therapy.



CPAP MACHINE – RESMED

AirSense 10 Autoset

ResMed

Available colors:

- black

- Automatic
- Humidifier (optional)
- Very quiet

TECHNICAL PARAMETERS	
dimensions	160 × 255 × 150 mm
noise level	26 dB
pressure range	4-20 mBar
weight	1.25 kg
RAMP function	AUTO mode or adjustable in the range 0-45 min
AC power supply	AC 100/240V, 50/60 Hz

Intended use:

An automatic device for the treatment of sleep apnea. Particularly recommended for people with a sensitive respiratory system. It has a comfort function and an AUTO-Ramp function which control the gentle way of air supply. The applied pre-heating function prevents the cooling of the respiratory tract during the initial operation of the device.



CPAP MACHINE – RESMED

AirSense 10 Autoset for Her



Available colors:

○ white

- Automatic
- Humidifier (optional)
- Very quiet

TECHNICAL PARAMETERS	
dimensions	160 × 255 × 150 mm
noise level	26 dB
pressure range	4-20 mBar
weight	1.25 kg
RAMP function	AUTO mode or adjustable in the range 0-45 min
AC power supply	AC 100/240V, 50/60 Hz

**Intended use:**

An automatic machine, in which, apart from all other functions characteristic for the devices of this company, a special algorithm adapted to female breathing was also used. Its use guarantees greater comfort thanks to slower and smaller fluctuations in the air pressure. The machine eliminates obstructive apnea, working in the algorithm adapted to the way women breathe up to 12 cm H₂O, and continuously controls and reacts appropriately to the existing air flow restrictions in the universal algorithm from 12 cm H₂O to 20 cm H₂O.

CPAP MACHINE – RESMED

AirSense 10 Elite



Available colors:

● black

- Constant airway pressure
- Humidifier (optional)
- Very quiet

TECHNICAL PARAMETERS	
dimensions	160 × 255 × 150 mm
noise level	26 dB
pressure range	4-20 mBar
weight	1.25 kg
RAMP function	AUTO mode or adjustable in the range 0-45 min
AC power supply	AC 100/240V, 50/60 Hz

**Intended use:**

A constant airway pressure machine for the treatment of sleep apnea. Particularly recommended for people with a sensitive respiratory system. It has a comfort function and an AUTO-Ramp function which control a gentle way of air supply. The applied pre-heating function prevents the cooling down of the respiratory tract during the initial operation of the device.

CPAP MACHINE – RESMED

AirMini



Available colors:

☐ white

- Automatic
- Very small
- Very quiet
- Controlled only via smartphone

TECHNICAL PARAMETERS	
dimensions	136 × 84 × 52 mm
noise level	29 dB
pressure range	4-20 mBar
weight	0.3 kg
RAMP function	AUTO mode or adjustable in the range 0-45 min
AC power supply	AC 100/240V, 50/60 Hz

Intended use:

An automatic machine intended for people who travel a lot. Its small size and weight make it convenient to move, and at the same time the therapy can be used in different places without losing the comfort of a peaceful sleep.

NOTE: Only the following mask models can be used: AirFit N20, AirFit F20, AirFit F30, and the P10 AirMini starter kit and the N30 AirMini starter kit.



BIPAP MACHINE – YUWELL

YH-830



Available colors:

☐ white

- Integrated humidifier
- Quiet
- Possibility to set different operating modes

TECHNICAL PARAMETERS	
dimensions	272.5 × 181 × 93 mm
noise level	30 dB
pressure: ipap	4-30 mBar
pressure: epap	4-21 mBar
weight	1.5 kg
RAMP function	Adjustable in the range 0-45 min
AC power supply	AC 100/240V, 50/60 Hz

Intended use:

Intended for the treatment of central apnea. This machine is distinguished by a modern design, light weight and easy operation. It indicates all the necessary parameters in the treatment of sleep apnea, including the AHI index, run time, mask leaks, which can be read directly on the monitor of the device. In addition, it is possible to save the therapy data by using an SD card and read a detailed report using a computer program.



MASK – RESMED

AirFit N30



Available colors:

● gray

- Subnasal
- Silicone

TECHNICAL PARAMETERS	
sizes	S and M
intended use	For all CPAP and BiPAP machines

Intended use:

Intended for CPAP machine users who can breathe freely through the nose and do not sleep at night with their mouth open. Due to the fact that the air is administered directly into the nostrils, it is intended rather for people with therapeutic pressure below 10 mBar.



MASK – RESMED

AirFit N20



Available colors:

● gray

- Nasal
- Silicone
- Magnetic headgear clips

TECHNICAL PARAMETERS	
sizes	S, M and L
intended use	For all CPAP and BiPAP machines

Intended use:

Intended for CPAP machine users who can breathe freely through the nose and do not sleep at night with their mouth open.



MASK – RESMED

AirFit N30i



Available colors:

● gray

- Subnasal
- Silicone
- Top-of-the-head air tube connection

TECHNICAL PARAMETERS	
sizes	S and M
intended use	For all CPAP and BiPAP machines

Intended use:

Intended for CPAP machine users who can breathe freely through the nose and do not sleep at night with their mouth open. Due to the fact that the air is administered directly into the nostrils, it is intended rather for people with therapeutic pressure below 10 mBar.



MASK – RESMED

AirFit F20



Available colors:

● gray

- Facial (oronasal)
- Silicone
- Magnetic clips

TECHNICAL PARAMETERS	
sizes	S, M and L
intended use	For all CPAP and BiPAP machines

Intended use:

Intended for CPAP machine users who can't breathe freely through the nose or open their mouth during sleep. It covers the nose and mouth, ensuring proper air supply.



MASK – RESMED

AirFit F30




- Facial (oronasal)
- Silicone
- Magnetic clips

TECHNICAL PARAMETERS	
sizes	S and M
intended use	For all CPAP and BiPAP machines

Intended use:

Intended for CPAP machine users who can't breathe freely through the nose or open their mouth during sleep. It covers the nose and mouth, ensuring proper air supply.



Available colors:

● gray

MASK – RESMED

AirFit F30i



- Facial (oronasal)
- Silicone
- Magnetic clips
- Top-of-the-head air tube connection

TECHNICAL PARAMETERS	
sizes	S and M
intended use	For all CPAP and BiPAP machines

Intended use:

Intended for CPAP machine users who can't breathe freely through the nose or open their mouth during sleep. The top of the mask covers only the lower part of the nose, the remaining part covers the mouth, ensuring proper air supply.



Available colors:

● gray

MASK – RESMED

AirFit P10



- Intranasal
- Silicone
- Size pillows included

TECHNICAL PARAMETERS	
sizes	S, M and L
intended use	For all CPAP and BiPAP machines

Intended use:

Intended for CPAP machine users who can breathe freely through the nose and do not sleep at night with their mouth open. Due to the fact that the air is administered directly into the nostrils, it is intended rather for people with therapeutic pressure below 10 mBar.



Available colors:

● blue

● pink

MASK – YUWELL YN-03

yuwell

Available colors:

● gray

- Nasal
- Silicone
- Hook fasteners

TECHNICAL PARAMETERS	
sizes	S, M and L
intended use	For all CPAP and BiPAP machines

Intended use:

Intended for CPAP machine users who can breathe freely through the nose and do not sleep at night with their mouth open.



MASK – YUWELL YF-02

yuwell

Available colors:

● gray

- Facial (oronasal)
- Silicone
- Hook fasteners

TECHNICAL PARAMETERS	
sizes	S, M and L
intended use	For all CPAP and BiPAP machines

Intended use:

Intended for CPAP machine users who can't breathe freely through the nose or open their mouth during sleep. It covers the nose and mouth, ensuring proper air supply.







**ORTHOPEDIC
PRODUCTS**



To say that we are only **effective and inexpensive** is like saying nothing about us.

The OMAX brand is:



aesthetic and elegant
at the same time



well-arranged
but not boring



simple in its form
but efficient



joyful, energetic and,
above all, active

OMAX are orthopedic products of the highest quality, providing maximum protection and comfort of use.

OMAX answers the question: What products would we like to use? Simple, practical, comfortable and effective. Thanks to this, now you can take advantage of our offer of orthopedic products.

What can we offer you?

We know that the needs of our customers are diverse, that's why we have created 3 product lines based on different requirements:



is a group of orthopedic products called "soft orthopedics", intended for everyday protection of joints during strenuous activities and ensuring slight stabilization during post-traumatic recovery. Made mostly of warming premium grade neoprene.



is a range of professional orthopedic equipment meeting the requirements of reimbursement by the National Health Fund. This line includes braces used in the process of perioperative and postoperative treatment, as well as in chronic and advanced diseases. These are products that allow the user to feel safe and comfortable. They combine top quality materials and technological solutions. They ensure the highest level of protection and comfort of use.



is a line of products designed for athletes and active people. Some of these products meet the requirements of reimbursement by the National Health Fund. These particular products are intended for the protection of muscles and joints during exercise and practicing various sports, as well as for post-traumatic treatment. They combine excellent quality with an aesthetic design. The products in this range are light, airy and flexible. They fit the body perfectly, providing excellent stabilization and protection.

Who are we?

As the REHA FUND company, the owner of the OMAX brand, we have been operating on the Polish and European medical market since 1993. We are a manufacturer and distributor of medical devices. We listen to the needs of our customers and adapt our products to them. You can find our products in medical shops all over Poland.

Find out more on www.omaxortopedia.pl

Foot and ankle



O-1180 Boot Brace "AIR WALKER" long

O-1181 Boot Brace "AIR WALKER" short



- Hard shell
- Soft, removable inner liner
- Two-chamber air system
- Integrated pump and vent valve
- Removable front protection panel
- Anti-slip, cradle sole
- Adjustable and removable toe protector
- Universal model for the right and left leg
- Available in a long and short version

Size	Shoe size
S	35-38
M	38-40
L	40-43
XL	43-46

Intended use:

Limb stabilization and immobilization in the case of ankle sprains and fractures, stabilization of the ankle joint, stress fractures and fractures of the metatarsal bones.



O-1180



O-1181

O-1182 Boot Brace

"Air Walker Basic" long



O-1183 Boot Brace

"Air Walker Basic" short

- Lightweight splint design
- Two-chamber air system
- Integrated pump and vent valve
- Soft, removable inner liner
- Anti-slip, cradle sole
- Universal model for the right and left leg
- Available in a long and short version

Size	Shoe size
S	35-38 cm
M	38-40 cm
L	40-43 cm
XL	43-46 cm

O-1182



Intended use:

Limb stabilization and immobilization in the case of ankle sprains and fractures, stabilization of the ankle joint, stress fractures and fractures of the metatarsal bones.

O-1183



O-1140 Ankle Stabilizer

"Achimax"



- Highly flexible compression knitted fabric
- Anatomically shaped silicone pad along the Achilles tendon
- Built-in heel pad
- Additional heel riser included
- Available in right and left leg versions

Size	Ankle circumference
S	24-26 cm
M	26-28 cm
L	28-31 cm
XL	31-34 cm
XXL	34-37 cm



Intended use:

Prevention and support in the treatment of injuries within the Achilles tendon.

O-1120 Neoprene Ankle Stabilizer

"Ligamax Ankle"



- Premium grade neoprene
- Open design
- System of cross straps with hook and loop fastening
- Universal model for the right and left leg

Size	Ankle circumference above the ankle
1	up to 32 cm
2	up to 38 cm



Intended use:

Used as a preventive measure in sport and when returning to physical activity after an injury.

O-1160 Ankle Stabilizer

"Malleomax"



- Rigid plastic side reinforcements
- Comfortable and breathable fabric
- Anatomically shaped, lace-up construction
- Open heel
- Elastic tightening straps
- Universal model for the right and left leg

Size	Ankle circumference
S	15-22 cm
M	22-30 cm
L	30-36 cm

Intended use:

Ankle joint stabilization and immobilization in the case of instability, mild to moderate sprains and dislocations of the joint, post-traumatic recovery.



O-1161 Ankle Stabilizer

"Malleomax Pro"



- Rigid plastic side reinforcements
- Comfortable and breathable fabric
- Anatomically shaped design
- Rotary adjustable compression system
- Open heel
- Elastic tightening straps
- Universal model for the right and left leg

Size	Ankle circumference
S	15-22 cm
M	22-30 cm
L	30-36 cm

Intended use:

Ankle joint stabilization and immobilization in the case of instability, mild to moderate sprains and dislocations of the joint, post-traumatic recovery.



Knee

O-1270 Postoperative Knee Brace "T-Max Knee"



- Telescopic aluminum side splints
- Fastening straps with quick-release buckles
- Removable anti-slip foam cushions
- Lightweight
- Wide range of length adjustment: 46-67 cm
- Extension angle adjustment: 0°-90° (the angle changes every 10°)
- Flexion angle adjustment: 0°-120° (the angle changes every 10°)
- Universal model for the right and left leg
- Universal size

Intended use:

Postoperative immobilization, ligament damage, joint protection in the process of post-traumatic rehabilitation.



O-1240 Long Closed Knee Brace
"Stabilmax Short"



- Compression and warming fabric
- Semi-open design
- Built-in peri-patellar silicone ring
- Flexion and extension angle adjustment (0°, 15°, 30°, 60°, 90°)
- Brace length 40 cm
- Universal model for the right and left leg

Size	Knee circumference
S	30–34cm
M	34.5–38cm
L	38.5–42cm
XL	42.5–46cm
XXL	46.5–50cm
3XL	50.5–54cm
4XL	54.5–58cm

Intended use:

Support and lateral-medial stabilization of the knee joint, protection in the process of post-traumatic rehabilitation, postoperative treatment and in disease entities requiring control of the motion range in the joint.



O-1241 Long Open Knee Brace
"Stabilmax Open Short"



- Compression and warming fabric
- Open design with back hook and loop fastening
- Built-in peri-patellar silicone ring
- Flexion and extension angle adjustment (0°, 15°, 30°, 60°, 90°)
- Brace length 40 cm
- Universal model for the right and left leg

Intended use:

Support and lateral-medial stabilization of the knee joint, protection in the process of post-traumatic rehabilitation, postoperative treatment and in disease entities requiring control of the motion range in the joint.

Size	Knee circumference
S	30–34cm
M	34.5–38cm
L	38.5–42cm
XL	42.5–46cm
XXL	46.5–50cm
3XL	50.5–54cm
4XL	54.5–58cm



O-1242 Short Closed Knee Brace

"Stabilmax Long"



- Compression and warming fabric
- Semi-open design
- Built-in peri-patellar silicone ring
- Flexion and extension angle adjustment (0°, 15°, 30°, 60°, 90°)
- Brace length 30 cm
- Universal model for the right and left leg

Intended use:

Support and lateral-medial stabilization of the knee joint, protection in the process of post-traumatic rehabilitation, postoperative treatment and in disease entities requiring control of the motion range in the joint.

Size	Knee circumference
S	30–34cm
M	34.5–38cm
L	38.5–42cm
XL	42.5–46cm
XXL	46.5–50cm
3XL	50.5–54cm
4XL	54.5–58cm



O-1243 Short Open Knee Brace

"Stabilmax Open Long"



- Compression and warming fabric
- Open design with back hook and loop fastening
- Built-in peri-patellar silicone ring
- Flexion and extension angle adjustment (0°, 15°, 30°, 60°, 90°)
- Brace length 30 cm
- Universal model for the right and left leg

Intended use:

Support and lateral-medial stabilization of the knee joint, protection in the process of post-traumatic rehabilitation, postoperative treatment and in disease entities requiring control of the motion range in the joint.

Size	Knee circumference
S	30–34cm
M	34.5–38cm
L	38.5–42cm
XL	42.5–46cm
XXL	46.5–50cm
3XL	50.5–54cm
4XL	54.5–58cm



O-1210 Flexible Knee Support with Hinges "Genumax Plus"

- Highly elastic compression knitted fabric
- Silicone ring to control the kneecap
- Side stays with hinges
- Two straps with hook and loop fastening
- Anatomical shape with comfort zones
- Anti-slip silicone inner lining
- Universal model for the right and left leg

Intended use:

Protection and stabilization of the joint in the case of mild to moderate instability, degenerative disease, inflammation.

Size	Knee circumference at the level of the kneecap
S	34-36 cm
M	36-38 cm
L	38-41 cm
XL	41-44 cm
XXL	44-47 cm



O-1211 Flexible Knee Support with Silicone Ring and Stays "Genumax"

- Highly elastic compression knitted fabric
- Silicone ring to control the kneecap
- Side underwires
- Anatomical shape with comfort zones
- Anti-slip silicone inner lining
- Universal model for the right and left leg

Intended use:

Joint protection and stabilization in the case of mild joint instability, chronic and post-traumatic inflammation.

Size	Knee circumference at the level of the kneecap
S	34-36 cm
M	36-38 cm
L	38-41 cm
XL	41-44 cm
XXL	44-47 cm



O-1201 Neoprene Knee Stabilizer

"Ligamax Knee"



- Premium grade neoprene
- Open design with hook and loop fastening
- Built-in silicone ring
- Universal model for the right and left leg

Size	Knee circumference at the level of the kneecap
1	32-43 cm
2	44-54 cm



Intended use:

Used as a preventive measure in sport and when returning to physical activity after an injury.

O-1212 Infra-patellar Stabilizer

"Rotulamax"



- Premium grade neoprene
- Open design with hook and loop fastening
- Anatomically shaped
- Cylindrical insert supporting the patella
- Universal model for the right and left leg
- Universal size



Intended use:

Support and stabilization in the case of patellar injuries, in the initial phase of instability of the knee joint.

Abdomen



O-1400 Abdominal Binder "24cm" "Ventrimax"



- Elastic compression fabric
- Open design with hook and loop fastening
- Binder height 24 cm

Available colors:

- O-1400 beige
- O-1400/CZA black

Intended use:

Protection and support after abdominal surgery, cesarean section, prevention of postoperative hernia.

Size	Waist circumference
S	81-90 cm
M	91-100 cm
L	101-110 cm
XL	111-120 cm
XXL	121-130 cm



O-1401 Breathable Abdominal Binder "24cm" "Ventrimax Air"



- Elastic perforated compression fabric
- Open design with hook and loop fastening
- Binder height 24 cm

Available colors:

- O-1401 beige
- O-1401/CZA black

Intended use:

Protection and support after abdominal surgery, cesarean section, prevention of postoperative hernia.

Size	Waist circumference
S	81-90 cm
M	91-100 cm
L	101-110 cm
XL	111-120 cm
XXL	121-130 cm



O-1404 Ostomy Abdominal Binder "Ventrimax Stoma"



- Elastic compression fabric
- Open design with hook and loop fastening
- Binder height 24 cm

Available colors:

- O-1404 beige
- O-1404/CZA black

Intended use:

Protection and support in the case of intestinal fistula, hernia in the area of the fistula, in the rehabilitation of patients with a stoma in the postoperative period.

Size	Waist circumference
S	81-90 cm
M	91-100 cm
L	101-110 cm
XL	111-120 cm
XXL	121-130 cm



O-1408 Single-sided Inguinal-hernia Belt "Axilmax"



- Elastic compression fabric
- Wide hip belt
- Open design with hook and loop fastening
- Inguinal strap with hook and loop fastening
- One inguinal pad
- Available in a right (O-1408P) and left (O-1408L) version

Available colors:

- O-1408 beige
- O-1408/CZA black

Intended use:

Support of the inguinal hernia, preventively after surgery, in the case of a weakened myofascial system.

Size	Waist circumference
S	88-97 cm
M	98-107 cm
L	108-117 cm
XL	117-125 cm



O-1409 Double-sided Inguinal-hernia Belt "Axilmax Duo"



- Elastic compression fabric
- Wide hip belt
- Open design with hook and loop fastening
- Two inguinal straps with hook and loop fastening
- Two inguinal pads

Available colors:

- O-1409 beige
- O-1409/CZA black

Intended use:

Support of the inguinal hernia, preventively after surgery, in the case of a weakened myofascial system.

Size	Waist circumference
S	88-97 cm
M	98-107 cm
L	108-117 cm
XL	117-125 cm



Spine



O-1461 High Semi-rigid Corset with Arm Supports "Lumbomax Spina"



- Elastic compression fabric
- System of airy channels
- Two-piece design
- Main belt with hook and loop fastening, two removable stays and two underwires
- Hook and loop fastened panel with arm supports and two removable stays
- Cross straps
- Centrally fastened tightening straps
- Main belt height 32 cm + back panel

Available colors:

- O-1461 beige
- O-1461/CZA black



Intended use:

Support for the spine in the case of degenerative changes, inflammation, overloading pain of the thoracolumbosacral spine.



Size	Waist circumference
S	70-80 cm
M	80-90 cm
L	90-100 cm
XL	100-110 cm
XXL	110-120 cm
3XL	120-135 cm
4XL	135-150 cm

O-1460 High Semi-rigid Corset with Tightening Straps

"Lumbomax Stabil"



- Elastic compression fabric
- System of airy channels
- Hook and loop fastening
- Two removable stays
- Two underwires
- Cross straps
- Centrally fastened tightening straps
- Corset height 32 cm



Available colors:

- O-1460 beige
- O-1460/CZA black

Intended use:

Support for the spine in the case of degenerative changes, inflammation, overloading and pain syndrome of the thoracolumbosacral spine.


Size	Waist circumference
S	70-80 cm
M	80-90 cm
L	90-100 cm
XL	100-110 cm
XXL	110-120 cm
3XL	120-135 cm
4XL	135-150 cm

O-1431 Lumbo-sacral Corset with Tightening Straps

"Lumbomax Plus"



- Elastic compression fabric
- System of airy channels
- Hook and loop fastening
- Two removable stays
- Two underwires
- Cross straps
- Centrally fastened tightening straps
- Corset height 26 cm



Available colors:

- O-1431 beige
- O-1431/CZA black

Intended use:

Support for the spine in the case of degenerative changes, inflammation, overloading and pain syndrome of the lumbosacral spine.

Size	Waist circumference
S	70-80 cm
M	80-90 cm
L	90-100 cm
XL	100-110 cm
XXL	110-120 cm
3XL	120-135 cm
4XL	135-150 cm

O-1430 Lumbo-sacral Corset

"Lumbomax"



- Elastic compression fabric
- System of airy channels
- Hook and loop fastening
- Four underwires
- Cross straps
- Corset height 26 cm



Available colors:

- O-1430 beige
- O-1430/CZA black

Intended use:

Support for the spine in the case of degenerative changes, inflammation, overloading and pain syndrome of the lumbosacral spine.

Size	Waist circumference
S	70-80 cm
M	80-90 cm
L	90-100 cm
XL	100-110 cm
XXL	110-120 cm
3XL	120-135 cm
4XL	135-150 cm

O-1432 Lumbo-sacral Corset with Silicone Pad

"Lumbomax Skin"



- Highly elastic compression knitted fabric
- Hook and loop fastening
- Two underwires
- Hook and loop fastened massaging pad
- Corset height 26 cm



Intended use:

Support for the spine in the case of degenerative changes, inflammation, overloading pain of the lumbosacral spine.

Size	Waist circumference
S	70-80 cm
M	80-90 cm
L	90-100 cm
XL	100-110 cm
XXL	110-120 cm

O-1433 Lumbo-sacral Corset with Tightening Straps
"Lumbomax Disc"

- Highly elastic compression knitted fabric
- Hook and loop fastening
- Four underwires
- Hook and loop fastened back panel
- Centrally fastened tightening straps
- Corset height 26 cm


Size	Waist circumference
S	70-80 cm
M	80-90 cm
L	90-100 cm
XL	100-110 cm
XXL	110-120 cm

Intended use:

Support for the spine in the case of degenerative changes, inflammation, overloading pain of the lumbosacral spine.



Shoulder



O-1500 Arm Sling



- Durable waterproof fabric
- Adjustable shoulder strap
- Anatomically shaped
- Soft shoulder pad
- Buckle closure
- Thumb holder
- Universal model for the right and left arm

Intended use:

Support and protection in the case of upper limb paresis, conservative or postoperative treatment.



Size	Length of the forearm measured in flexion (from elbow to wrist)
S	18-24 cm
M	21-27 cm
L	28-33 cm
XL	34-40 cm

O-1501 Mesh Arm Sling



- Soft and breathable mesh fabric
- Adjustable shoulder strap
- Anatomically shaped
- Soft shoulder pad
- Buckle closure
- Thumb holder
- Universal model for the right and left arm

Available colors:

- O-1501/N blue
- O-1501/R pink
- O-1501/CZA black

Intended use:

Support and protection in the case of upper limb paresis, conservative or postoperative treatment.



Size	Length of the forearm measured in flexion (from elbow to wrist)
XS	13-18 cm
S	18-24 cm
M	21-27 cm
L	28-33 cm
XL	34-40 cm

- Made of breathable perforated polyurethane foam
- Open design
- Anatomically shaped
- Hook and loop fastened
- Double-sided insertion function
- Universal model for the right and left arm


Size	Waist circumference
S	65-75 cm
M	75-90 cm
L	90-105 cm
XL	105-120 cm
XXL	120-135 cm
3XL	135-145 cm



Intended use:

Support and stabilization in the case of shoulder joint dislocations, in acute pain syndromes of the shoulder, in the treatment of injuries of the shoulder joint, shoulder girdle and elbow joint, as well as in postoperative stabilization of the shoulder joint.

Elbow



O-1670 Range of Motion Elbow Brace "T-Max Elbow"



- Single-axis telescopic splint
- Soft lining cushions
- Adjustable straps with hook and loop fastening
- Lightweight
- Length adjustment range: 33-43 cm
- Extension angle adjustment: -15° to 90° (the angle changes every 15°)
- Flexion angle adjustment: 0° - 105° (the angle changes every 15°)
- Universal size
- Available in a right/left version

Intended use:

Postoperative immobilization, ligament damage, joint protection in the process of post-traumatic rehabilitation.



O-1620 Elbow Support with Silicone Pad
"Epimax"



- Highly elastic compression knitted fabric
- Anatomically shaped silicone massage pads
- Anatomical shape with comfort zones
- Additional strap with hook and loop fastening
- Universal model for the right and left arm

Size	Circumference of the elbow joint
S	22-25 cm
M	25-28 cm
L	25-31 cm
XL	31-34 cm
XXL	34-37 cm

Intended use:

Prevention and support in the treatment of elbow injuries, elbow pain syndromes.



O-1601 Neoprene Elbow Support
"Ligamax Elbow"




- Premium grade neoprene
- Open design with hook and loop fastening
- Universal model for the right and left arm
- Universal size

Intended use:

Used as a preventive measure in sport and when returning to physical activity after an injury.



Wris and Hand



O-1742 Elastic Wrist and Thumb Support "Manimax T"



- Highly elastic compression knitted fabric
- Adjustable straps with hook and loop fastening
- Anatomical shape with comfort zones
- Anatomically shaped palm splint
- Additional underwires stiffening the back of the hand
- Flexible aluminum thumb support splint
- Wide wristband
- Available in a right/left version

Intended use:

Joint stabilization and protection in the case of: injuries of the radiocarpal joint, overloading, degenerative and inflammatory conditions.



Size	Wrist circumference
S	15-17 cm
M	17-19 cm
L	19-21 cm
XL	21-23 cm
XXL	23-25 cm

O-1743 Elastic Wrist Support "Manimax"



- Highly elastic compression knitted fabric
- Open design
- Adjustable straps with hook and loop fastening
- Built-in comfort zones
- Anatomically shaped palm splint
- Additional underwires stiffening the back of the hand
- Wide wristband
- Available in a right/left version

Intended use:

Joint stabilization and protection in the case of: injuries of the radiocarpal joint, overloading, degenerative and inflammatory conditions.



Size	Wrist circumference
S	15-17 cm
M	17-19 cm
L	19-21 cm
XL	21-23 cm
XXL	23-25 cm

O-1740 Forearm and Hand Brace with Thumb
 "Manumax T"



- Made of cotton fabric and soft foam
- Anatomically shaped, removable hand splint
- Additional stays stiffening the back of the hand
- Flexible aluminum thumb support splint
- Adjustable straps with hook and loop fastening
- Available in a right/left version



Intended use:

Joint stabilization and protection in the case of injuries of the radiocarpal joint, overloading, degenerative and inflammatory conditions.

Size	Wrist circumference
S	14-18 cm
M	19-22 cm
L	23-28 cm

O-1741 Forearm and Hand Brace
 "Manumax"



- Made of cotton fabric and soft foam
- Open design
- Anatomically shaped, removable hand splint
- Additional stays stiffening the back of the hand
- Adjustable straps with hook and loop fastening
- Universal model for the right and left hand



Intended use:

Joint stabilization and protection in the case of injuries of the radiocarpal joint, overloading, degenerative and inflammatory conditions.

Size	Wrist circumference
S	14-16 cm
M	17-19 cm
L	20-23 cm

O-1760 Short Wrist Brace

"Manumax Uni Short"



- Lightweight open design
- Anatomically shaped aluminum stay
- Additional universal thumb stay
- Three tightening straps to ensure proper fit to the hand
- Brace length 20 cm
- Universal model for the right and left hand
- Universal size



Intended use:

Joint stabilization and protection in the case of wrist and thumb injuries, wrist and thumb instability, overloading, degenerative and inflammatory conditions.

O-1761 Long Wrist Brace

"Manumax Uni Long"



- Lightweight open design
- Anatomically shaped aluminum stay
- Additional universal thumb stay
- Three tightening straps to ensure proper fit to the hand
- Brace length 25 cm
- Universal model for the right and left hand
- Universal size



Intended use:

Joint stabilization and protection in the case of wrist and thumb injuries, wrist and thumb instability, overloading, degenerative and inflammatory conditions.

O-1770 Multi-finger Brace
"Digimax"



- Made of velour fabric and neoprene
- Glove design
- Hook and loop fastening
- Four channels for stiffening splints
- One finger splint included
- Available in a right/left version



Intended use:

Stabilization in the case of finger extensor injuries, some types of fractures of the basal and middle phalanges.

Size	Wrist circumference
S	17-19 cm
M	19.5-22 cm
L	22.5-25 cm

O-1771 Finger Brace
"Digimax Plus"



- Made of breathable fabric
- Aluminum stiffening splint
- Adjustable hook and loop fastening around the wrist
- Additional stabilizing straps with hook and loop fastening
- Universal model for the right and left hand



Intended use:

Joint stabilization and protection in the case of radiocarpal joint injuries, overloading, degenerative and inflammatory conditions.

Size	Wrist circumference
S	13-17 cm
M	17.5-21 cm

O-1720 Neoprene Wrist Support
"Ligamax Wrist"



- Premium grade neoprene
- Open short design with hook and loop fastening
- Universal model for the right and left hand
- Universal size

Intended use:

Used as a preventive measure in sport and when returning to physical activity after an injury.



O-1721 Neoprene Wrist Support with Thumb Loop
"Ligamax Palm"



- Premium grade neoprene
- Open design with hook and loop fastening
- Universal model for the right and left hand
- Universal size

Intended use:

Used as a preventive measure in sport and when returning to physical activity after an injury.



Neck

O-1800 Schanza Type Orthopedic Collar



- Foam covered with cotton fabric
- Anatomically shaped
- Open design with hook and loop fastening
- Additional washable cover
- Available in three heights: 8 cm, 10 cm, 12 cm

Intended use:

Used in the case of degenerative and post-traumatic changes in the cervical spine, pain and overloading syndromes of the neck and shoulders.



Size	Neck circumference
S	34-37 cm
M	37-42 cm
L	42-47 cm
XL	47-54 cm

O-1840 Philadelphia Type Tracheostomy rigid collar



- Made of polyethylene foam
- Anatomically shaped
- PVC reinforcement
- Two-piece design
- Hook and loop fastening
- Equipped with a tracheostomy opening

Intended use:

Used in the case of pain and radicular syndromes of the cervical spine, injuries of the cervical spine, in post-accident immobilization.



Size	Neck circumference
S	34-36 cm
M	36-38 cm
L	38-40 cm





**REHABILITATION
PRODUCTS**

ABS REHABILITATION BALL

PILO-1045/-1055/-1065/-1075

- Durable break-resistant PVC material
- Maximum safe load up to 300 kg
- Available sizes: 45 cm, 55 cm, 65 cm, 75 cm
- Pump included

reference number	size	height	color
PILO-1045	45 cm	140 – 154 cm	yellow
PILO-1055	55 cm	155 – 169 cm	red
PILO-1065	65 cm	170 – 184 cm	blue
PILO-1075	75 cm	above 185 cm	black



Intended use:

ABS ball for rehabilitation and general development exercises. Perfect for exercises that correct posture defects, improve fitness, strengthen muscles and increase joint mobility. Used for balancing and stretching exercises, relieving back pain and excessive muscle tension.

MASSAGE REHABILITATION BALL

PILO-2045/-2055/-2065/-2075

- Durable break-resistant PVC material
- Maximum safe load up to 300 kg
- Available sizes: 45 cm, 55 cm, 65 cm, 75 cm
- Pump included

reference number	size	height	color
PILO-2045	45 cm	140 – 154 cm	yellow
PILO-2055	55 cm	155 – 169 cm	red
PILO-2065	65 cm	170 – 184 cm	blue
PILO-2075	75 cm	above 185 cm	black



Intended use:


ABS ball for rehabilitation and general development exercises. Perfect for exercises that correct posture defects, improve fitness, strengthen muscles and increase joint mobility. Used for balancing and stretching exercises, relieving back pain and excessive muscle tension. Equipped with soft protrusions on the entire surface, which additionally stimulate skin receptors and provide micro-massage and muscle relaxation.

REHABILITATION BALL "MINI"

PILO-1020/-1025

- Durable break-resistant PVC material
- Available sizes: 20 cm, 25 cm
- Pump included

reference number	size	color
PILO-1020	20 cm	pink
PILO-1025	25 cm	blue



Intended use:

"MINI" ball for rehabilitation and general development exercises. Perfect for isometric exercises to prevent muscle atrophy, improve fitness and increase joint mobility. It is ideal for the rehabilitation of posture defects and spinal disorders. Used for pilates and yoga.


"PEANUT" ABS REHABILITATION BALL

PILO-1090

- Durable break-resistant PVC material
- Maximum safe load up to 250 kg
- Size: 90 cm × 40 cm
- Color ●
- Pump included

Intended use:

"ABS peanut" ball for rehabilitation and general development exercises. Perfect for exercises that correct posture defects, improve fitness, strengthen muscles and increase joint mobility. Its specially profiled shape allows for movement only in two directions (forward and backward), which gives greater control of movement and more pronounced stability of the exercises.



REHABILITATION BALL WITH SPIKES


PILO-3006/-3007/-3008/-3009/-3010

- Made of the highest quality odorless and non-toxic PVC material
- Available sizes: 6 cm, 7 cm, 8 cm, 9 cm, 10 cm

reference number	size	color
PILO-3006	6 cm	●
PILO-3007	7 cm	●
PILO-3008	8 cm	●
PILO-3009	9 cm	●
PILO-3010	10 cm	●

Intended use:

Ball with spikes intended for rehabilitation, hand and foot exercises and a body massage. Equipped with protrusions on the entire surface, which additionally stimulate skin receptors, improve blood circulation and microcirculation. They reduce tension, relax muscles and de-stress.



REHABILITATION BALL WITH SPIKES SOFT

– SET


PILO-3011

- Made of the highest quality odorless and non-toxic PVC material
- The set includes 3 balls in sizes 7 cm, 8 cm and 9 cm

reference number	size	color
PILO-3011	7 cm	●
	8 cm	●
	9 cm	●

Intended use:

Ball with spikes intended for rehabilitation, hand and foot exercises and a body massage. Equipped with protrusions on the entire surface, which additionally stimulate skin receptors, improve blood circulation and microcirculation. They reduce tension, relax muscles and de-stress.



"CUCUMBER" REHABILITATION BALL

PILO-3015

- Made of the highest quality odorless and non-toxic PVC material
- Available size: 15 cm x 7 cm

reference number	color
3015/1	green
3015/2	black
3015/3	red

Intended use:

"Cucumber" ball with spikes intended for rehabilitation, hand and foot exercises and a body massage. Equipped with protrusions on the entire surface, which additionally stimulate skin receptors, improve blood circulation and microcirculation. They reduce tension, relax muscles and de-stress. Its specially profiled shape ensures a better fit to the body.



"PEANUT" REHABILITATION BALL


PILO-3018

- Made of the highest quality odorless and non-toxic PVC material
- Available size: 17 cm x 8 cm

reference number	color
PILO-3018/1	blue
PILO-3018/2	black
PILO-3018/3	red

Intended use:

"Peanut" ball with spikes intended for rehabilitation, hand and foot exercises and a body massage. Equipped with protrusions on the entire surface, which additionally stimulate skin receptors, improve blood circulation and microcirculation. They reduce tension, relax muscles and de-stress. Its specially profiled shape ensures a better fit to the body.



"HEDGEHOG" SENSORY HEMISPHERE

PILO-2010

- Durable break-resistant PVC material
- Size: 16.5 cm diameter; 8 cm height

reference number	color
PILO-2010/1	green
PILO-2010/2	black
PILO-2010/3	red

Intended use:

Hemisphere for foot rehabilitation, motor coordination and balance exercises. Equipped with protrusions on the raised surface that stimulate deep sensation, improve blood circulation and microcirculation. They reduce tension, relax muscles and de-stress.



"JELLY" REHABILITATION BALL

PILO-5015/-5020/-5030

- Made of the highest quality thermoplastic rubber
- Size: 5 cm diameter
- Available 3 degrees of elasticity

reference number	degrees of elasticity	color
PILO-5015	soft	yellow
PILO-5020	medium	orange
PILO-5030	hard	blue

Intended use:

"Jelly" rehabilitation ball intended for exercises strengthening hand and forearm muscles. It improves dexterity and motor skills. Its use has an anti-stress and calming effect.



"EGG" REHABILITATION BALL

PILO-6015/-6020/-6030

- Made of the highest quality thermoplastic rubber
- Size: 5.8 cm x 4.4 cm
- Available 3 degrees of elasticity

reference number	degrees of elasticity	color
PILO-6015	soft	yellow
PILO-6020	medium	orange
PILO-6030	hard	blue

Intended use:

"Egg" rehabilitation ball intended for exercises strengthening hand and forearm muscles. It improves dexterity and motor skills. Its use has an anti-stress and calming effect.



REHABILITATION RING WITH SPIKES

PILO-3030/-3040/-3050

- Made of the highest quality silicone rubber
- Size: 7 cm outer diameter
- Available degrees of elasticity: green (soft), blue (medium), orange (hard)

reference number	color
PILO-3030	green
PILO-3040	blue
PILO-3050	orange

Intended use:



Rehabilitation ring with spikes intended for exercises strengthening hand and forearm muscles. It improves dexterity, motor skills and improves hand gripping functions.



SENSORY DISK

PILO-2034

- Durable break-resistant PVC material
- Size: 34 cm diameter
- Pump included

reference number	color
PILO-2034/1	
PILO-2034/2	





Intended use:

Sensory disc intended for post-traumatic rehabilitation, balance, sensorimotor and corrective exercises. Perfect for exercises that correct posture defects, improve fitness, strengthen muscles and increase joint mobility. Equipped with protrusions on the upper surface, which additionally stimulate skin receptors, provide micromassage and muscle relaxation. While sitting, it relieves the lumbar spine and forces the correct position.

"BALANCE" SENSORY DISK

PILO-2038

- Durable break-resistant PVC material
- Size: 38 cm diameter
- Pump included

reference number	color
PILO-2038/1	
PILO-2038/2	



Intended use:

"BALANCE" sensory disc intended for post-traumatic rehabilitation, balance, sensorimotor and corrective exercises. Perfect for exercises that correct posture defects, improve fitness, strengthen muscles and increase joint mobility. Equipped with protrusions on the upper surface, which additionally stimulate skin receptors, provide micromassage and muscle relaxation. While sitting, it relieves the lumbar spine and forces the correct position.

"PILOTAPE" SILK KINESIOLOGY TAPE

PILO-9010

- backing material made of light artificial silk
- width: 5 cm
- length: 5 m
- waterproof and quick drying
- hypoallergenic acrylic adhesive coating
- latex-free

reference number	color
9010/1	orange
9010/2	black
9010/3	pink
9010/4	blue
9010/5	green

Intended use:

Used to support the treatment and rehabilitation of injured body parts. Used in post-traumatic conditions of muscles and joints to relieve swelling, as well as to correct posture defects. Properly applied, it helps relieve pain, reduce muscle tension and support the work of the circulatory system. The thickness and stretch of the tape is similar to skin thickness and elasticity, allowing for a full range of motion.



"PILOTAPE KIDS" KINESIOLOGY TAPE

PILO-9020

- backing material made of light artificial silk
- width: 5 cm
- length: 5 m
- waterproof and quick drying
- hypoallergenic acrylic adhesive coating
- latex-free

reference number	color
PILO-9020/1	flamingo and pineapple pattern
PILO-9020/2	fox and forest pattern

Intended use:

Used to support the treatment and rehabilitation of injured body parts. Used in post-traumatic conditions of muscles and joints to relieve swelling, as well as to correct posture defects. Properly applied, it helps relieve pain, reduce muscle tension and support the work of the circulatory system. The thickness and stretch of the tape is similar to skin thickness and elasticity, allowing for a full range of motion.





AUXILIARY AND SANITARY AIDS

Manufacturer of sanitary and self-care aids

In our product range you can find all the necessary accessories that will certainly prove invaluable during the routine of the people struggling with disabilities. We produce not only **bed pans** or **urinal bottles**, but also **toilet and bathtub seats, reachers, shower stools** and **chairs** or **shampoo basins**.

Even simple pieces of equipment, such as a shower chair or toilet safety rail, provide great support on the way to self-reliance.

Our products meet very high quality standards, so you can be sure that they will prove to be reliable and fully functional!

Eating, washing yourself or using the toilet does not have to be a challenge. Choosing the Reha Fund offer, you can quickly find out that everyday life with a disability does not have to be burdensome. This equipment not only increases the level of independence of people with disabilities, but also allows them to enjoy the comfort of life.

We are fully aware how difficult it is to choose the right solutions, so our highly qualified staff will willingly help the customer in purchasing the ideal self-care aids, tailored to the individual needs of the future user.

TOILET SEAT – BRUNO 801

RF-801

Available colors:

● gray

- serves as a portable toilet and toilet extension
- made of steel
- foldable frame
- anti-slip rubber tips

TECHNICAL PARAMETERS	
width	55 cm
total height	72~82 cm
total length	72 cm
seat height	42-52 cm
seat depth	65 cm
backrest height	33 cm
max. safe load	100 kg
weight	6.6 kg

Intended use:

The toilet seat functions as a portable toilet. The seats can be used at the patient's bedside. It helps people with reduced mobility and the elderly in performing basic physiological activities. It also makes it easier to care for them. The barriers make it easier to get up from the seat and sit on it.



TOILET WHEELCHAIR – MARCO 802

RF-802

Available colors:

● black

- serves as a portable toilet
- can be used in the shower
- made of aluminum
- flip-up armrests
- swing away and detachable footrests
- seat and backrest made of waterproof material

TECHNICAL PARAMETERS	
total width	55 cm
seat width	46 cm
seat depth	40 cm
wheel diameter	13 cm
seat height	50.5 cm
total length	85 cm
height adjustment	95.5-98 cm
backrest height adjustment	45.5-48 cm
max. safe load	120 kg
weight	13.5 kg

Intended use:

The toilet wheelchair functions as a portable toilet. It ensures the maintenance of personal hygiene. It can be used at the patient's bedside or in the bathroom. It can also be used to transport the user by an assistant.



TOILET SEAT – MICHAEL 805

RF-805

- serves as a portable toilet
- self-assembly kit
- anti-slip rubber tips

Available colors:

● gray

TECHNICAL PARAMETERS	
total length	64 cm
total width	51 cm
total depth	59.5 cm
seat width	49.5 cm
seat depth	41 cm
seat height adjustment range	41-56.3 cm
max. safe load	120 kg
chair weight	6 kg

Intended use:

The toilet seat functions as a portable toilet. It can be used at the patient's bedside. It helps people with reduced mobility and the elderly in performing basic physiological activities. It also makes it easier to care for them. The barriers make it easier to get up from the seat and sit on it.



TOILET SEAT – LOUIS 840

RF-840

- 2 seats: comfortable – removable, with toilet function – with an opening, swing away and detachable
- height-adjustable legs with a safe lock
- anti-slip rubber tips

Available colors:

● black

● blue

TECHNICAL PARAMETERS	
max. safe load	130 kg
weight	8 kg
total width	56 cm
total height	80-95 cm
depth	52-57 cm
seat width	43 cm
seat depth	42 cm
seat height	40-55 cm
toilet hole size	24 x 21 cm

Intended use:

The toilet seat functions as a portable toilet. It can be used at the patient's bedside. It helps people with limited mobility (permanent or temporary) and the elderly in performing basic physiological activities. The seat can also be used as a toilet top and an ordinary chair to sit on. The seat facilitates patient care. The barriers make it easier to get up from the seat and sit on it. It is intended for users with limited mobility, but with sufficient seating ability. The seat is useful in cases of: patients with limited control of the trunk, people with balance disorders, dizziness, people with impaired function of the lower limbs, patients with significantly limited mobility of the lumbar spine, hip and knee joints, paralysis.



PERSONAL HYGIENE CHAIR – OLIWIER 810

RF-810

- facilitates the care of intimate areas
- with a U-shaped cut-out
- suitable for use in the shower
- adjustable height
- with a backrest for ensuring stability

TECHNICAL PARAMETERS	
total width	48.5 cm
seat width	48.5 cm
total height	72.5-82.5 cm
total length	47 cm
seat height	39-49 cm
seat depth	39 cm
backrest height	36 cm
max. safe load	110 kg
weight	2.6 kg

Intended use:

The chair allows people with limited mobility, people after operations and seniors to maintain daily personal hygiene, while staying significantly independent. A special cut-out also facilitates taking care of intimate hygiene. Thanks to this equipment, the risk of falling or other injuries that are possible in the shower is maximally reduced.

Available colors:

○ white



PERSONAL HYGIENE STOOL – WILSON 811

RF-811

- facilitates the care of intimate areas
- with a U-shaped cut-out
- suitable for use in the shower
- adjustable height

TECHNICAL PARAMETERS	
total width	48.5 cm
seat width	48.5 cm
total height	41-51 cm
total length	41 cm
seat height	39-49 cm
seat depth	33.5 cm
max. safe load	110 kg
weight	2 kg

Intended use:

The stool allows people with limited mobility, people after operations and seniors to maintain daily personal hygiene, while staying significantly independent. A special cut-out also facilitates taking care of intimate hygiene. Thanks to this equipment, the risk of falling or other injuries that are possible in the shower is maximally reduced.

Available colors:

○ white



SHOWER CHAIR – NATHAN 820

RF-820

- backrest for stability
- with side openings for easy getting up and sitting
- suitable for use in the shower
- adjustable height

TECHNICAL PARAMETERS	
total width	50 cm
seat width	50 cm
total height	75-85 cm
total length	48 cm
seat height	39-49 cm
seat depth	35 cm
backrest height	37 cm
max. safe load	110 kg
weight	2.9 kg

Intended use:

The chair allows people with limited mobility, people after operations and seniors to maintain daily personal hygiene, while staying significantly independent. Thanks to this equipment, the risk of falling or other injuries that are possible in the shower is maximally reduced.

Available colors:

○ white



SHOWER STOOL – KEVIN 821

RF-821

Available colors:

☐ white

- with side openings for easy getting up and sitting
- suitable for use in the shower
- adjustable height

TECHNICAL PARAMETERS	
total width	50 cm
seat width	50 cm
total height	41-51 cm
total length	39 cm
seat depth	31 cm
seat height	39-49 cm
max. safe load	110 kg
weight	1.9 kg



Intended use:

The stool allows people with limited mobility, people after operations and seniors to maintain daily personal hygiene, while staying significantly independent. Thanks to this equipment, the risk of falling or other injuries that are possible in the shower is maximally reduced.

ROUND SHOWER STOOL – CAMERON 822

RF-822

Available colors:

☒ blue

☐ white

- universal shower stool
- compact round form, ideal even for a smaller shower

TECHNICAL PARAMETERS	
max. safe load	136 kg
stool weight	2.5 kg
seat height adjustment	40-57.5 cm
seat diameter	32 cm

Intended use:

The stool allows people with limited mobility, people after operations and seniors to maintain daily personal hygiene, while staying significantly independent. Thanks to this equipment, the risk of falling or other injuries that are possible in the shower is maximally reduced.



2-IN-1 BATHTUB STOOL AND STEP – JOHN 823

RF-823

Available colors:

○ white

- compact size
- suitable for use inside or outside a bathtub

TECHNICAL PARAMETERS	
max. safe load	120 kg
weight	1.5 kg
seat height adjustment	20.2-25.2 cm
seat width	41.6 cm
seat depth	22.6 cm
leg spacing lengthwise	47.8 cm
leg spacing crosswise	31.6 cm
leg angle lengthwise	107 degrees
leg angle crosswise	100 degrees

Intended use:

The stool/step allows people with limited mobility to maintain daily personal hygiene, while staying significantly independent. Thanks to this equipment, the risk of falling or other injuries that are possible in the bathtub is significantly reduced. The stool/step makes it easier to use the bathtub, helps in daily hygiene, and increases safety. The product makes it easier for people with limited mobility to get in and out of the bathtub. It can also be used as a stable support during bathing in the bathtub.



ROTATING ROUND SHOWER STOOL – RYAN 824

RF-824

Available colors:

● blue

- soft, swivel seat
- tray under the seat for storing cosmetics
- adjustable height

TECHNICAL PARAMETERS	
max. safe load	150 kg
stool weight	2.6 kg
seat height adjustment	42-54 cm
seat diameter	35 cm

Intended use:

The stool allows people with limited mobility to maintain daily personal hygiene, while staying significantly independent. Thanks to this equipment, the risk of falling or other injuries that are possible in the shower is significantly reduced.



BATHTUB SEAT – EMMA 851

RF-851

- adjustable bathtub fixation system
- made of easy-to-wash PVC material
- universal and easy to install

Available colors:

☐ white

TECHNICAL PARAMETERS	
max. safe load	181 kg
weight	2.14 kg
mounting width adjustment	500-600 mm
seat width	68.5 cm
seat depth	22.5 cm
total height	11.5 cm
slat height	1.8 cm
slat width	5 cm

Intended use:

This product provides a stable seat while bathing in a standard bathtub. It assists in getting in and out of the tub, in sitting down and/or getting up in the tub. It is intended for the elderly and people with limited mobility (temporary or permanent) caused by e.g. stroke, hemorrhage, traffic accident, bone fragility, injuries of limbs, hips, surgery.



ADJUSTABLE BATHTUB SEAT – AVA 853

RF-853

- deep seat that allows to sit down lower in the bathtub
- adjustable width

Available colors:

☐ white

TECHNICAL PARAMETERS	
max. safe load	150 kg
total width	71-81 cm
mounting spacing	28.3 cm
seat width	41.6 cm
seat length	22.6 cm
height of the lowered part	14.5 cm
weight	2 kg

Intended use:

This product provides a stable seat when bathing in a standard bathtub. It assists in getting in and out of the tub, in sitting down and/or getting up in the tub. It is intended for the elderly and people with limited mobility (temporary or permanent) caused by e.g. stroke, hemorrhage, traffic accident, bone fragility, injuries of limbs, hips, surgery.



SHOWER SEAT – SOPHIA 854

RF-854

- stable steel frame
- wall-mounted
- mounting kit attached
- foldable legs
- adjustable height

Available colors:

☐ white

TECHNICAL PARAMETERS	
max. safe load	113 kg
height adjustment	40-55 cm
seat width	49.5 cm
seat depth	30 cm
weight	4.8 kg

Intended use:

The product provides a stable seat while showering. It assists in everyday care activities and sitting up while in the shower. It is intended for the elderly and people with limited mobility (temporary or permanent) caused by e.g. stroke, hemorrhage, traffic accident, bone fragility, injuries of limbs, hips, surgery.



BATHTUB SEAT – OLIVIA 855

RF-855

Available colors:

○ white

- made of easy-to-wash plastic
- universal and easy to install
- ergonomic handle and opening for cosmetics

TECHNICAL PARAMETERS	
max. safe load	100 kg
length	71 cm
seat width	32 cm
seat height	10 cm
weight	2.5 kg
total height	18.5 cm

Intended use:

The product provides a stable seat when bathing in a standard bathtub. It assists in getting in and out of the tub, in sitting down and/or getting up in the tub. It is intended for the elderly and people with limited mobility (temporary or permanent) caused by e.g. stroke, hemorrhage, traffic accident, bone fragility, injuries of limbs, hips, surgery.



BATHTUB SEAT – DIANA 856

RF-856

Available colors:

○ white

- made of easy-to-wash plastic
- universal and easy to install
- ergonomic handle and opening for cosmetics
- with a seat

TECHNICAL PARAMETERS	
max. safe load	110 kg
seat width	35 cm
mounting spacing	46-68 cm
weight	2.7 kg
total height	9.5 cm
total length	76 cm

Intended use:

The product provides a stable seat when bathing in a standard bathtub. It assists in getting in and out of the tub, in sitting down and/or getting up in the tub. It is intended for the elderly and people with limited mobility (temporary or permanent) caused by e.g. stroke, hemorrhage, traffic accident, bone fragility, injuries of limbs, hips, surgery.



RAISED TOILET SEAT – SARA 870

RF-870

- increases the height of a toilet
- easy to install
- with a lid
- anatomical shape
- fits all standard toilets
- tool-free modern fixing system with build-in screws
- height 10 or 15 cm

TECHNICAL PARAMETERS		
MODELS		
reference number	RF-870/10	RF-870/15
max. safe load	100 kg	100 kg
height	100 mm	150 mm
length	390 mm	390 mm
seat width	360 mm	360 mm
weight	1.4 kg	1.6 kg

Intended use:

The raised toilet seat is used to increase the height of the toilet to the most comfortable level for the user. This makes it easier to sit down and stand up. It is intended for all people who find it difficult to move around and use the toilet in a traditional way, e.g. for patients after hip surgery or injury, for seniors.

Available colors:

☐ white



RAISED TOILET SEAT – KATE 871

RF-871

- increases the height of a toilet
- easy to install with locking clips
- with a lid
- anatomical shape
- height 11 or 15 cm

TECHNICAL PARAMETERS		
MODELS		
reference number	RF-871/11	RF-871/15
max. safe load	150 kg	150 kg
front height	117 mm	147 mm
back height	147 mm	177 mm
back height with clips	207 mm	237 mm
seat width	367 mm	367 mm
seat depth	401 mm	401 mm
opening width	208 mm	208 mm
opening depth	264 mm	264 mm
weight	1.05 kg	1.32 kg

Intended use:

The raised toilet seat is used to increase the height of the toilet to the level most comfortable for the user. This makes it easier to sit down and stand up. It is intended for all people who find it difficult to move around and use the toilet in a traditional way, e.g. for patients after hip surgery or injury, for seniors.

Available colors:

☐ white



KIDNEY DISH

RF-894

- easy to keep clean
- suitable for cold sterilization and for sterilization in dry air at 130°C

TECHNICAL PARAMETERS		
models	Small	Big
capacity	300 ml	700 ml
dimensions	20.5 x 11.4 x 3.85 cm	28.2 x 14.65 x 4.5 cm

Intended use:

The product is intended for use in medical units and at home to service bedridden people, in order to alleviate the effects of injuries and the course of diseases. The product is used for washing injuries, post-traumatic care and other hygienic procedures. When the product is used by elderly or immobilized people, third party assistance is required.

Available colors:

☐ white



DRINKING CUP

RF-895

Available colors:

○ white

- comes with 2 interchangeable lids – for drinking fluids and for semi-liquid foods
- equipped with 2 easy grip handles
- suitable for cold sterilization (e.g. by using ethanol)

TECHNICAL PARAMETERS	
capacity	200 ml
dimensions	14.4 cm × dia. 7.1 cm

Intended use:

The product is intended for use in medical units and at home to service bedridden people, in order to alleviate the effects of injuries and the course of diseases. The drinking cup is equipped with two lids that enable serving meals to patients in a lying position. The lid with the small opening should be used for drinking liquids; the lid with the large opening should be used for semi-liquid food (mush). When the product is used by elderly or immobilized people, third party assistance is required.



SHAMPOO BASIN

RF-960

Available colors:

○ white

- air inflatable
- draining pipe fixed
- special cut-out for the neck
- intended for washing patient's hair in a lying position

TECHNICAL PARAMETERS	
length	53 cm
width	61 cm
height	20 cm
weight	0.83 kg
capacity	10 liters

Intended use:

Intended for bedridden or immobilized people.



SHAMPOO BASIN WITH WATER BAG

RF-961

Available colors:

○ white

- air inflatable
- special cut-out for the neck
- draining pipe fixed
- shower set included
- intended for washing patient's hair in a lying position

TECHNICAL PARAMETERS	
dimensions when inflated (length × width × height)	53 × 61 × 20 cm
dimensions when deflated and folded (length × width × height)	22 × 25 × 4 cm
drain pipe length	118 cm
PVC material thickness	0.35 cm
weight	1.6 kg

Intended use:

Intended for bedridden or immobilized people.



STRAIGHT BATHROOM GRAB BAR – ALEX 910

RF-910

- ideal for people who need extra support in everyday toilet activities
- universal use with all kinds of sanitary devices
- made of high quality plastic
- wall-mounted
- mounting kit included



Available colors:

○ white

TECHNICAL PARAMETERS				
MODELS				
reference number	RF-910/30	RF-910/40	RF-910/45	RF-910/60
max. safe load	100 kg			
length	300 mm	400 mm	450 mm	600 mm
tube diameter	350 mm			
distance from the wall	850 mm			
weight	0.22 kg	0.29 kg	0.36 kg	0.43 kg

Intended use:

The bathroom grab bar serves as a support element. It can be used in a damp room.

STRAIGHT BATHROOM GRAB BAR – CEZAR 911

RF-911

- ideal for people who need extra support in everyday toilet activities
- universal use with all kinds of sanitary devices
- durable and stable steel construction
- powder coated
- mounting kit included



Available colors:

○ white

TECHNICAL PARAMETERS					
model	RF-911/20	RF-911/30	RF-911/40	RF-911/50	RF-911/60
bar length	20 cm	30 cm	40 cm	50 cm	60 cm
tube diameter	3.2 cm				
distance from the wall	11 cm				
max. safe load	135 kg				
weight	0.35 kg	0.40 kg	0.55 kg	0.7 kg	0.85 kg

Intended use:

The grab bar helps people with limited mobility in maintaining balance and moving, increases the safety of people using sanitary facilities, is suitable for installation in the bathroom. It is intended for use in apartments, nursing homes and facilities occupied by disabled people. The product is mainly used to prevent falls.

BATHTUB SUPPORT RAIL – DENNIS 900

RF-900

- made of high quality steel
- powder coated
- strong and durable
- tool-free and screw-free assemble
- universal for various types of bathtubs

TECHNICAL PARAMETERS	
max. safe load	180 kg
weight	2.82 kg
total height	48.5 cm
width adjustment	8.2~17.2 cm


Intended use:

The support rail helps in maintaining balance and moving for people with limited mobility, increases the safety of people using sanitary facilities, suitable for mounting in the bathroom on the bathtub without using tools. It is intended for use in apartments, nursing homes and facilities occupied by disabled people.



Available colors:

○ white



FLIP-UP TOILET SUPPORT RAIL – DANIEL 920

Available colors:


○ white

- durable steel construction
- installed mainly beside the toilet
- powder coated
- mounting kit included

TECHNICAL PARAMETERS	
max. safe load	136 kg
weight	2.85 kg
total length	75 cm
total height	23 cm

Intended use:

The support rail helps maintain balance and move for people with reduced mobility, increases the safety of people using sanitary facilities, is suitable for mounting at the washbasin, toilet or another desired place in the bathroom. It is intended for use in apartments, nursing homes, facilities occupied by disabled people.



TOILET SAFETY FRAME – JUDY 940

Available colors:

● gray

- aluminum structure, fixed to the toilet seat lid mounting bolts
- anti-slip rubber tips
- height and width adjustable

TECHNICAL PARAMETERS	
max. safe load	113 kg
weight	1.55 kg
height adjustment	64-74 cm
width adjustment	53-63 cm

Intended use:

The safety frame is intended to help when getting up and sitting on the toilet for the elderly, the disabled and people having problems with maintaining balance.



MALE URINAL BOTTLE

RF-890

- equipped with a grip and a secure cup
- suitable for cold sterilization, for dry air sterilization at 130°C or steam sterilization in an autoclave at 126°C

TECHNICAL PARAMETERS	
capacity	1000 ml
dimensions	28 × 12 × 11.3 cm

Intended use:

The product is intended for use in medical units and at home to service bedridden people, in order to alleviate the effects of injuries and the course of diseases. The product is used to collect urine. When the product is used by elderly or immobilized people, third party assistance is required.

Available colors:

☐ white



FEMALE URINAL BOTTLE

RF-891

- has a specially profiled opening
- suitable for cold sterilization, for dry air sterilization at 130°C or steam sterilization in an autoclave at 126°C

TECHNICAL PARAMETERS	
capacity	700 ml
dimensions	31 × 12 × 11.5 cm

Intended use:

The product is intended for use in medical units and at home to service bedridden people, in order to alleviate the effects of injuries and the course of diseases. The product is used to collect urine. When the product is used by elderly or immobilized people, third party assistance is required.

Available colors:

☐ white



SLIPPER BED PAN

RF-892

- equipped with a comfortable handle and a secure fitting lid to prevent spillage
- suitable for cold sterilization, for dry air sterilization at 130°C or steam sterilization in an autoclave at 126°C

TECHNICAL PARAMETERS	
capacity	3000 ml
dimensions	12.05 × 29.15 × 54.75 cm

Intended use:

The product is intended for use in medical units and at home to service bedridden people in order to alleviate the effects of injuries and the course of diseases. The product is compatible with paper pulp inserts and can also be used as their base. The use of the product in combination with disposable cartridges made of cellulose pulp significantly reduces the possibility of infection spreading. The product is used to collect feces and urine in bedridden people. When the product is used by elderly and/or immobilized people, third party assistance is required.

Available colors:

☐ white



FLAT SLIPPER BED PAN

RF-896

- flat design with a handle
- possibility of completing with disposable cellulose inserts
- suitable for dry air sterilization at 130°C or steam sterilization in an autoclave at 126°C

TECHNICAL PARAMETERS	
capacity	2000 ml
dimensions	8 × 27 × 48 cm

Intended use:

The product is intended for use in medical units and at home to service bedridden people in order to alleviate the effects of injuries and the course of diseases. The product is compatible with paper pulp inserts and can also be used as their base. The use of the product in combination with disposable cartridges made of cellulose pulp significantly reduces the possibility of infection spreading. The product is used to collect feces and urine in bedridden people. When the product is used by elderly and/or immobilized people, third party assistance is required.



Available colors:

- white
- green
- blue

DISPOSABLE SLIPPER BED PAN

MED-02

- disposable cellulose slipper bed pan
- can be used as an insert for the RF-896 slipper bed pan

TECHNICAL PARAMETERS	
dimensions	36.5 × 27 × 8.5 cm (at the deepest point)
capacity	2000 ml

Intended use:

The disposable slipper bed pan is intended for people who cannot go to the toilet on their own. The product makes it possible to satisfy one's physiological needs in a lying or semi-sitting position. It is used in medical units and at home to service bedridden people. When used by the elderly or immobilized people, third party assistance is required.



Available colors:

- gray

DAILY URINE COLLECTION CONTAINER – TULIP

RF-893

- equipped with a spill stopper lid
- suitable for cold sterilization, for dry air sterilization at 130°C or steam sterilization in an autoclave at 126°C

TECHNICAL PARAMETERS	
dimensions	19 cm × Ø18.5 cm
capacity	2500 ml

Intended use:

The product is intended for use in health care units for daily urine collection. When the product is used by elderly or immobilized people, third party assistance is required.



Available colors:

- white

UPPER LIMB REHABILITATION DEVICE 881

RF-881

- helps restore the motion range of the upper body
- simple and reliable design
- simple assembly without the need to interfere with the infrastructure (over the door, without screws)

TECHNICAL PARAMETERS	
weight	0.75 kg

Intended use:

The device is intended for exercises of the upper limbs, shoulders and upper body. It helps in rebuilding the range of motion and increasing the motion/mobility range in the joints of the upper limbs and upper torso. The resistance for one limb is the weight of the other limb.



Available colors:

● gray

BED BACKREST 950

RF-950

- adjustable almost to a sitting position (max. 80°)
- 5-stage adjustment

TECHNICAL PARAMETERS	
max. safe load	75 kg
angle adjustment range	0°-80°
dimensions	47 × 63 cm
weight	2.64 kg

Intended use:

The bed backrest is intended for people who require a higher torso position while lying.



Available colors:

● blue

DENTURE CUP

DNT

- box with a convenient lid
- strainer for easy drainage
- non-medical device

TECHNICAL PARAMETERS	
dimensions	9 × 8 × 6 cm

Intended use:

Reusable denture box for use in nursing, rehabilitation, elderly care and personal use at home.



Available colors:

● green

SHORT REACHER

RF-750

- for picking up small items without bending down
- ended with a magnet to attract objects e.g. a coin or jewelry

TECHNICAL PARAMETERS	
reference number	RF-750
model name	Short reacher
length	66 cm
max. load	1 kg

Intended use:

The reacher is intended for people who have problems with hand mobility, most often caused by spinal diseases, accidents, which are followed by difficulty in bending. The reacher helps pick up small objects without having to bend down or to shorten the distance. The jaws with a built-in magnet also allow one to pick up small metal objects such as pins, needles, paper clips, etc.



Available colors:

● gray

LONG REACHER

RF-751

- for picking up small items without bending down
- ended with a magnet to attract objects e.g. a coin or jewelry

TECHNICAL PARAMETERS	
reference number	RF-751
model name	Long reacher
length	76 cm
max. load	1 kg

Intended use:

The reacher is intended for people who have problems with hand mobility, most often caused by spinal diseases, accidents, which are followed by difficulty in bending. The reacher helps pick up small objects without having to bend down or to shorten the distance. The jaws with a built-in magnet also allow one to pick up small metal objects such as pins, needles, paper clips, etc.



Available colors:

● gray

ORTHODONTIC BOX

702_R

- box with a convenient lid
- wide range of colors
- non-medical device

TECHNICAL PARAMETERS	
dimensions	9 × 8 × 6 cm

Intended use:

The orthodontic box is intended to protect against bacteria or mechanical damage to orthodontic appliances and removable retention appliances.



Available colors:

- blue
- green
- yellow
- pink
- orange
- red





PILL ORGANIZERS AND ACCESSORIES



Manufacturer of pill organizers and accessories

We have produced a line of pill organizers and accessories for patients who take multiple doses of medication throughout the day, for people who regularly supplement their diet or for those who travel a lot.

The division into times of day and week is an ideal daily help in organizing medicines. Thanks to the use of handy organizers, the patient is able to take medications regularly and remember about the appropriate times for administering the next dose.

DAILY PILL ORGANIZER (TWICE-A-DAY)

RF-350

- enables the division of medications into two daily portions
- made of polypropylene
- non-medical device

TECHNICAL PARAMETERS	
length	5.3 cm
depth	3.9 cm
height	2 cm

Intended use:

Used to divide the doses of medications.



DAILY PILL ORGANIZER (THREE-A-DAY)

RF-351

- enables the division of medications into three daily portions
- made of polypropylene
- non-medical device

TECHNICAL PARAMETERS	
length	7.8 cm
depth	3.9 cm
height	2 cm

Intended use:

Used to divide the doses of medications.



DAILY PILL ORGANIZER (FOUR-A-DAY)

RF-352

- enables the division of medications into four daily portions
- made of polypropylene
- non-medical device

TECHNICAL PARAMETERS	
length	10.3 cm
depth	3.9 cm
height	2 cm

Intended use:

Used to divide the doses of medications.



WEEKLY PILL ORGANIZER FOR DOSING ORAL PILLS


RF-353

- enables the division of medications into 4 doses for a 7 day period
- made of polypropylene
- non-medical device

TECHNICAL PARAMETERS	
length	14 cm
depth	9.4 cm
height	2.8 cm

Intended use:

Used to divide the doses of medications.



WEEKLY PILL ORGANIZER (ONCE-A-DAY)

RF-354

- enables the division of medications into daily doses for a 7 day period
- made of polypropylene
- non-medical device

TECHNICAL PARAMETERS	
length	16.2 cm
depth	3.9 cm
height	2 cm

Intended use:

Used to divide the doses of medications.



WEEKLY PILL ORGANIZER (TWICE-A-DAY)

RF-355

- enables the division of medications into 2 daily doses for a 7 day period
- made of polypropylene
- non-medical device

TECHNICAL PARAMETERS	
length	16.3 cm
depth	8.2 cm
height	2 cm

Intended use:

Used to divide the doses of medications.



WEEKLY PILL ORGANIZER (THREE-A-DAY)

RF-356

- enables the division of medications into 3 doses for a 7 day period
- made of polystyrene
- non-medical device

TECHNICAL PARAMETERS	
length	8.6 cm
depth	4.7 cm
height	14.3 cm

Intended use:

Used to divide the doses of medications.



WEEKLY PILL ORGANIZER (FOUR-A-DAY)

RF-357

- enables the division of medications into 4 doses for a 7 day period
- made of polystyrene
- non-medical device

TECHNICAL PARAMETERS	
length	11.2 cm
depth	4.7 cm
height	14.3 cm

Intended use:

Used to divide the doses of medications.



PILL CUTTER

RF-358

- intended to cut pills
- perfect for proper portioning
- made of polystyrene
- non-medical device

TECHNICAL PARAMETERS	
length	7.3 cm
depth	2.9 cm
height	2.6 cm
maximum tablet diameter	1.9 cm

Intended use:

It cuts pills, facilitates their precise portioning. Used also for storing medications.



PILL CRUSHER

RF-359

- crushes pills into powder to help with ingestion for children and adults
- comes with a container for storing crushed pills
- made of polystyrene
- non-medical device

TECHNICAL PARAMETERS	
length	6 cm
depth	3.6 cm
height	4.5 cm
opening diameter	2.3 cm

Intended use:

It crushes pills, making it easier for children and adults to take them. It has a container for storing medications.



TOOTHPASTE & CREAM TUBE SQUEEZER

RF-360

- helps precisely squeeze out the tube content
- made of polystyrene
- non-medical device

TECHNICAL PARAMETERS	
length	10.2 cm
depth	4.2 cm
height	3.6 cm
maximum width of tubes placed in the squeezer	6 cm

Intended use:

It enables accurate squeezing of the contents of the tube and improves the display of products. It can be used in the cosmetics, pharmaceutical and food industries.



TONGUE CLEANER

RF-361

- removes bacterial deposits and food residues from the tongue
- leaves fresh breath
- quick and easy to use
- if you have trouble flossing, use the floss holder
- made of polystyrene
- non-medical device

TECHNICAL PARAMETERS	
length	14.5 cm
depth	4 cm
height	0.7 cm

Intended use:

The cleaner removes bacterial deposits and food particles from the tongue.







ANTI-BEDSORE MATTRESSES

Manufacturer of anti-bedsore foam and pneumatic mattresses

Basing on the knowledge and experience gained over several decades on healthy sleep as well as prevention and treatment of ulcers, we have created a unique offer of anti-bedsore and therapeutic mattresses focused on the individual needs of the users. In our offer, the patient will find both bedsore mattresses, designed to reduce the pressure on the body as much as possible, while supporting the correct position, as well as anti-bedsore, therapeutic mattresses focused on supporting the correct body position while maximizing sleep comfort.

A foam or a pneumatic mattress?

This is a frequently asked question that does not have a clear answer. It all depends on the needs of the user. Foam mattresses provide effective anti-bedsore prevention and a stable base, while being less susceptible to mechanical damage. They often do not need an additional base and provide very good air circulation. On the other hand, alternating pressure mattresses provide definitely better anti-bedsore prevention because by changing the pressure in the chambers during one cycle, the patient's entire body is relieved. The pneumatic anti-bedsore mattress will be most effective for people at high risk of developing bedsores or already affected by them.

Stages of bedsores – an important parameter when choosing an anti-bedsore mattress

Stage I: Intact skin with non-blanchable redness of a localized area usually over a bony prominence. Darkly pigmented skin may not have visible blanching; its color may differ from the surrounding area. The site may be painful, hard, soft, warmer or cooler than the surrounding tissue.

Stage II: Partial thickness loss of dermis presenting as a shallow open ulcer with a red pink wound bed, without slough separating from healthy tissue. May also present as an intact or open/ruptured serum-filled blister. Presents as a shiny or dry shallow ulcer without slough or bruising.

Stage III: Full thickness tissue loss. Subcutaneous fat may be visible but a bone, tendon or muscle is not exposed. Slough may be present but does not obscure the depth of the tissue loss. It may include undermining and tunneling. The depth of a stage III bedsore varies by anatomical location.*


* European Pressure Ulcer Advisory Panel and National Pressure Ulcer Advisory Panel. Prevention and treatment of pressure ulcers: quick reference guide. Washington DC: National Pressure Ulcer Advisory Panel; 2009.

MEDICAL ALTERNATING PRESSURE PNEUMATIC MATTRESS WITH PUMP HF6001+HF62012



**Available
colors:**

● beige



- practical easy-to-use bubble mattress
- easy to install
- reduces pressure on blood vessels, adapts to the patient's body and position
- for people with an average or high risk of developing stage I bedsores
- for people who will spend more than 15 hours a day in bed
- for short-term care
- comes with a repair kit and patches
- made of highly elastic medical PVC material

MATTRESS HF6001

SPECIFICATION	
max. user weight	130 kg
mattress dimensions	200 cm × 90 cm × 7 cm
material	PVC (polyvinyl chloride) 0.3 mm
color	beige
mattress weight	1.8-2.1 kg
number of mattress cells	130
cycle time	12 min
warranty	18 months

PUMP HF62012

SPECIFICATION	
air output	5-7 l/min
pressure range	50-120 mmHg
power supply	220 V/50 Hz
capacity	5-7 l/min
fuses	1 A
power consumption	7 W
dimensions	23.5 (L) × 12 (H) × 9.5 (W) cm
warranty	18 months
pump weight	0.85 kg

Intended use:

A traditionally accepted method of bedsores prevention is to rotate the patient and to monitor risk areas regularly. The use of an alternating pressure mattress is a supplement to this treatment. At the same time, turning and observing the patient still needs to be considered a part of patient care, if needed.

PNEUMATIC MATTRESS

BIOFLOTE™ 2000

BioFlote™


Available colors:

● beige



- high-quality bubble mattress
- for people in need of treatment and prevention of bedsores up to stage I
- for people who are constantly in a lying position
- for short-term care
- for people in need of treatment and reduction of pain caused by prolonged lying
- large air bubbles that increase the anti-bedsores effect
- pump equipped with a precise pressure level indicator set according to the patient's weight
- pump equipped with hooks for mounting on the bed with anti-slip tips
- flaps to stabilize the patient
- repair kit including a patch with glue
- quiet pump – noise below 31 decibels
- flexible PVC hose

PUMP SPECIFICATION

power supply	230V / 50Hz
fuse	1 A
pump dimensions in cm	27 × 11 × 9,5
pump weight	1.2 kg
power consumption	7 W
cycle duration	10-12 min
pressure range	40-100 mmHg
pump capacity	5 liters/min

MATTRESS SPECIFICATION

type	Bubble (130 bubbles)
mattress dimensions in cm	200 × 90 × 6,5 – when inflated
mattress weight	1.9 kg
max. user weight	130 kg
material	Polyvinyl chloride (EN-71 compliant)
ambient temperature	Operating: 5°C~40°C Storage: 25°C~70°C Transport: – 25°C~70°C
ambient relative humidity	Operating: 15%~90% Storage: 10%~90% Transport: 10%~90%
atmospheric pressure	40-100 mmHg

Intended use:

This product is intended for those who: need treatment and prevention of stage I bedsores while optimizing patient comfort, permanently staying in one position, for short-term care at home; need treatment and reduction of pain caused by the pressure of the underneath surface on the soft tissue of the body with the recommendation of a doctor or physiotherapist. This product may be used by people who are qualified to perform general care procedures and who are adequately trained and know how to prevent and treat bedsores.

ANTI-BEDSORE MATTRESS

BIOFLOTE™ 3000

BioFlote™

Available colors:

● blue



- simple easy-to-use tubular mattress
- perfect for home conditions
- for people in need of treatment and prevention of stage II bedsores
- for people who are constantly in a lying position
- for people in need of treatment and reduction of pain caused by prolonged lying
- improves patient comfort
- built of 17 chambers
- spare chamber included
- pump equipped with a precise pressure level indicator set according to the patient's weight
- for long-term care
- pump equipped with hooks for mounting on the bed with anti-slip tips
- anchor openings to stabilize the mattress
- quiet pump – noise below 31 decibels
- single chamber stabilizing strap

PUMP SPECIFICATION

power supply	230V / 50Hz
protection	1 A
pump dimensions in cm	27 × 11 × 9.5
pump weight	1.2 kg
power consumption	7 W
cycle time	10~12 min
pressure range	40-100 mmHg
pump capacity	5 liters/min

MATTRESS SPECIFICATION

type	Tubular (17 air chambers)
mattress dimensions in cm	198 × 88 × 11 – when inflated
mattress weight	4.8 kg
max. user weight	130 kg
material	Nylon/polyvinyl chloride
ambient temperature	Operating: 5°C~40°C Storage: 25°C~70°C Transport: – 25°C~70°C
ambient relative humidity	Operating: 15%~90% Storage: 10%~90% Transport: 10%~90%
atmospheric pressure	40-100 mmHg

Intended use:

This product is intended for those who: need treatment and prevention of stage II bedsores while optimizing patient comfort, permanently staying in one position, for long-term home care and treatment of bedsores; need treatment and reduction of pain caused by the pressure of the surface on the soft tissue of the body with the recommendation of a doctor or physiotherapist. This product may be used by people who are qualified to perform general care procedures and who are adequately trained and know how to treat and prevent bedsores.


PNEUMATIC MATTRESS

BIOFLOTE™ 4000

BioFlote™

Available colors:

● blue



- advanced tubular mattress with additional features
- suitable for both home care and hospital use
- for people in need of treatment and prevention of stage II bedsores
- for people who are constantly in a lying position
- for people in need of treatment and reduction of pain caused by prolonged lying
- improves patient comfort
- built of 17 chambers
- spare chamber included
- pump equipped with a precise pressure level indicator set according to the patient's weight
- for long-term care
- pump equipped with transport quick plug
- pump equipped with anti-slip hooks
- CPR valve
- equipped with a cover fastened with press studs
- micro-ventilation system – ensures additional ventilation of the patient's skin
- pillow function – 3 chambers under the head can be inflated all the time
- double chamber stabilizing strap

PUMP SPECIFICATION

power supply	230V / 50Hz
fuse	1 A
pump dimensions in cm	28 × 15 × 9
pump weight	1.4 kg
power consumption	7 W
cycle time	12 min
pressure range	40-100 mmHg
pump capacity	8 liters/min

MATTRESS SPECIFICATION

type	Tubular (17 air chambers)
mattress dimensions in cm	190 × 88 × 11 – when inflated
mattress weight	4.85 kg
max. user weight	140 kg
material	Nylon/PVC
ambient temperature	Operating: 5°C~40°C Storage: – 25°C~70°C Transport: – 25°C~70°C
ambient relative humidity	Operating: 15%~90% Storage: 10%~90% Transport: 10%~90%
atmospheric pressure	40-100 mmHg

Intended use:

This product is intended for those who: need treatment and prevention of stage II bedsores while optimizing patient comfort, permanently staying in one position, for long-term home care and treatment of bedsores; need treatment and reduction of pain caused by the pressure of the surface on the soft tissue of the body with the recommendation of a doctor or physiotherapist. This product may be used by people who are qualified to perform general care procedures and who are adequately trained and know how to treat and prevent bedsores.


ANTI-BEDSORE MATTRESS

BIOFLOTE™ 5000

BioFlote™

Available colors:

● blue



- advanced tubular mattress with multiple features
- suitable for both home care and hospital use
- for people in need of treatment and prevention of stage III bedsores
- for people who are constantly in a lying position
- for people in need of treatment and reduction of pain caused by prolonged lying
- improves patient comfort
- pump equipped with a precise pressure level indicator set according to the patient's weight
- for long-term care
- built of 17 chambers
- spare chamber included
- pump equipped with quick-release transport plug
- pump equipped with anti-slip hooks
- CPR valve
- equipped with a zippered cover
- pillow function – 3 chambers under the head can be inflated all the time
- double chamber stabilizing strap
- static function inflates the entire mattress
- visual low pressure alarm informs caregivers about a possible mattress air leak

PUMP SPECIFICATION

power supply	230V / 50Hz
fuse	1 A
pump dimensions in cm	28 × 15 × 9
pump weight	2.2 kg
power consumption	7 W
cycle time	12 min
pressure range	30-60 mmHg
pump capacity	8 liters/min

MATTRESS SPECIFICATION

type	Tubular (17 air chambers)
mattress dimensions in cm	203 × 90 × 12.7 – when inflated
mattress weight	5.2 kg
max. user weight	160 kg
material	Nylon/PVC
ambient temperature	Operating: 5°C~40°C Storage: 25°C~70°C Transport: – 25°C~70°C
ambient relative humidity	Operating: 15%~90% Storage: 10%~90% Transport: 10%~90%
atmospheric pressure	30-60 mmHg

Intended use:

This product is intended for those who: need treatment and prevention of stage III bedsores while optimizing patient comfort, permanently staying in one position, for long-term home care and treatment of bedsores; need treatment and reduction of pain caused by the pressure of the surface on the soft tissue of the body with the recommendation of a doctor or physiotherapist. This product may be used by people who are qualified to perform general care procedures and who are adequately trained and know how to treat and prevent bedsores.

ANTI-BEDSORE MATTRESS

BIOFLOTE™ 50

BioFlote™



- intended to stabilize the pneumatic mattress against the underneath surface (e.g. bed)
- flat mattress profile
- foam density 28 kg/m³
- high foam hardness



TECHNICAL PARAMETERS

dimensions	200 × 90 × 5.5 cm
	200 × 85 × 5.5 cm
cover	without cover
	elastic waterproof

Intended use:

Anti-bedsore mattress with a flat profile, used as a base for a pneumatic mattress. Stabilizes the pneumatic mattress against the surface. Used to prevent the formation of stage I bedsores.

ANTI-BEDSORE MATTRESS

BIOFLOTE™ 200



BioFlote™



- standalone mattress
- mattress profile with a unique 3-zone shaped profile: convoluted ("egg crate") section in the head and feet area; longitudinal tube section in the pelvis area; the patient's head is higher than the legs
- foam density 25 kg/m³
- medium foam hardness

TECHNICAL PARAMETERS

dimensions	200 × 90 × 10–7 cm
	without cover
cover	elastic waterproof
	cotton-polyester
	waterproof



Intended use:

Anti-bedsore mattress with a unique 3-zone shaped profile, recommended for the prevention of stage I bedsores.

ANTI-BEDSORE MATTRESS


BIOFLOTE™ 250

BioFlote™



- standalone mattress
- convoluted ("egg crate") profile of the mattress facilitates air circulation and prevents bedsores
- foam density 25 kg/m³
- high foam hardness

TECHNICAL PARAMETERS	
dimensions	200 × 90 × 7 cm
	200 × 90 × 10 cm
cover	without cover
	elastic waterproof
	cotton-polyester
	waterproof



Intended use:
Anti-bedsores mattress with a convoluted ("egg crate") profile used to prevent the formation of stage I bedsores.

ANTI-BEDSORE MATTRESS


BIOFLOTE™ 300

BioFlote™



- standalone mattress
- made of two layers of foam
- the top layer of the mattress with a unique 3-zone shaped profile: convoluted ("egg crate") section in the head and feet area; longitudinal tube section in the pelvis area; the patient's head is higher than the legs
- the bottom layer of the foam stabilizes the mattress and reduces the pressure of the surface on the patient's body
- foam density 25 kg/m³
- medium foam hardness

TECHNICAL PARAMETERS	
dimensions	200 × 90 × 15–12 cm
	without cover
cover	elastic waterproof
	cotton-polyester
	waterproof



Intended use:
Anti-bedsores mattress with a unique 3-zone shaped profile, recommended for the prevention of stage I bedsores.


ANTI-BEDSORE MATTRESS

BIOFLOTE™ 350

BioFlote™



- standalone mattress
- made of two layers of foam
- the top layer of the mattress with a convoluted ("egg crate") profile
- the bottom layer of the foam stabilizes the mattress and reduces the pressure of the surface on the patient's body
- foam density – top layer: 25 kg/m³, bottom layer: 28 kg/m³
- high foam hardness



TECHNICAL PARAMETERS

dimensions	200 × 90 × 15 cm
	without cover
cover	elastic waterproof
	cotton-polyester
	waterproof

Intended use:

Anti-bedsore mattress with a convoluted ("egg crate") profile with an additional layer of stabilizing foam, used to prevent the formation of stage I bedsore.

ANTI-BEDSORE MATTRESS

BIOFLOTE™ 400 BASIC

BioFlote™



- standalone mattress
- flat mattress profile
- foam density 25 kg/m³
- low foam hardness



TECHNICAL PARAMETERS

dimensions	200 × 90 × 10 cm
	200 × 90 × 12 cm
	200 × 90 × 15 cm
cover	without cover
	elastic waterproof
	cotton-polyester
	waterproof

Intended use:

Anti-bedsore mattress with a flat profile, used to prevent the formation of stage I bedsore.

ANTI-BEDSORE MATTRESS

BIOFLOTE™ 600

BioFlote™



- standalone mattress
- flat mattress profile
- made of three layers of foam
- the top layer made of viscoelastic foam ensures use comfort
- the middle layer has transverse grooves to facilitate air circulation
- the bottom layer of the foam stabilizes the mattress and reduces the pressure of the surface on the patient's body
- foam density: top: 50 kg/m³, middle: 30 kg/m³, bottom: 28 kg/m³
- medium foam hardness

TECHNICAL PARAMETERS	
dimensions	200 × 90 × 15–12 cm
	without cover
cover	elastic waterproof
	cotton-polyester
	waterproof

Intended use:

Anti-bedsore mattress with a flat profile, used to prevent the formation of stage I bedsores.



ANTI-BEDSORE MATTRESS

BIOFLOTE™ 700

BioFlote™




- standalone or underlay mattress
- flat mattress profile
- foam density 18 kg/m³
- medium foam hardness

TECHNICAL PARAMETERS	
dimensions	200 × 90 × 7 cm
	200 × 90 × 10 cm
	without cover
cover	elastic waterproof
	cotton-polyester
	waterproof

Intended use:

Anti-bedsore mattress with a flat profile, used to prevent the formation of stage I bedsores.



ANTI-BEDSORE MATTRESS

BIOGOFER 15000

BioFlote™




- standalone mattress
- waffle-like profile of the mattress facilitating air circulation
- 152 independent segments
- foam density 30 kg/m³
- medium foam hardness

TECHNICAL PARAMETERS	
dimensions	200 × 90 × 10 cm
	200 × 90 × 12 cm
	200 × 90 × 15 cm
	without cover
cover	elastic waterproof
	cotton-polyester
	waterproof

Intended use:

Anti-bedsore mattress with a waffle-like profile, used to prevent the formation of stage I bedsores.



WATERPROOF ELASTIC COVER

15/P/EN

- waterproof
- reusable
- washable at 90°C
- suitable for foam and pneumatic mattresses
- elastic, with elastic band
- available color: white
- dimensions 200 × 90 × 15 cm

Intended use:

The elastic waterproof mattress cover with an elastic band is intended to protect the surface against contamination with physiological secretions.



COTTON-POLYESTER COVER

15/P/B

- reusable
- washable at 60°C
- suitable for foam mattresses
- dimensions: 200 × 90 × 7 cm, 200 × 90 × 10 cm, 200 × 90 × 12 cm, 200 × 90 × 15 cm

Intended use:

The cotton-polyester mattress cover is intended to protect the surface against contamination with physiological secretions.



WATERPROOF COVER

15/P/N

- waterproof
- reusable
- washable at 95°C
- suitable for foam mattresses
- dimensions: 200 × 90 × 7 cm, 200 × 90 × 10 cm, 200 × 90 × 12 cm, 200 × 90 × 15 cm

Intended use:

The waterproof vapor-permeable mattress cover is intended to protect the surface against contamination with physiological secretions.



REUSABLE HYGIENIC UNDERPAD


3045

Available colors:

- green
- blue

- waterproof
- inner material made of cotton-polyester terry knitted fabric
- outer waterproof material: 100% PVC
- reusable
- washable at 95°C
- resistant to disinfectants
- available colors: blue, green

TECHNICAL PARAMETERS	
dimensions	70 × 100 cm
	70 × 140 cm
	100 × 140 cm
	140 × 200 cm



Intended use:

This reusable hygienic underpad is a medical device intended to protect surfaces (e.g. mattress, bed) against contamination with physiological secretions (e.g. urine, blood). It protects against the occurrence of infections / infections as a result of prolonged contact of the patient with bodily secretions.

HYGIENIC UNDERPAD FOR HOME CARE

3046

Available colors:

- green
- blue
- white

- waterproof
- outer material made of cotton-polyester terry knitted fabric
- inner material – waterproof membrane: 100% polyurethane
- reusable
- washable at 90°C
- available colors: blue, green, white

TECHNICAL PARAMETERS	
dimensions	60 × 70 cm
	70 × 100 cm
	90 × 100 cm
	70 × 140 cm
	90 × 140 cm
	100 × 140 cm
	110 × 140 cm
	90 × 200 cm
	140 × 200 cm
	180 × 200 cm
	90 × 200 cm with elastic band
	140 × 200 cm with elastic band
	180 × 200 cm with elastic band



Intended use:

This hygienic underpad for home care is a medical device intended to protect surfaces (e.g. mattress, bed) against contamination with physiological secretions (e.g. urine, blood). It protects against the occurrence of infections resulting also from prolonged contact of the patient with bodily secretions.





ORTHOPEDIC PILLOWS FOR SLEEPING

Manufacturer of Morfeusz orthopedic pillows for sleeping

The acquired knowledge allowed us to create a wide range of **comfortable orthopedic pillows**, characterized by a carefully thought-out structure.

An incorrect sleeping position can lead to malformation of the cervical vertebrae and inappropriate muscle tension. As a consequence, it can cause pressure on the spinal cord, blood flow disorders, which is felt in the form of fatigue, neck and nape pain, numbness of the hands and dizziness. Therefore, a properly selected, **comfortable orthopedic pillow** adjusted to the body structure and sleeping position is crucial.

The foam that the orthopedic pillow is made of **has an open-chamber structure**. This allows free air circulation and eliminates the discomfort caused by head sweating.

Our orthopedic sleeping pillow does not cause allergies!

The MORFEUSZ orthopedic pillows have covers made of **100% cotton knitted fabric**. Cotton is a natural non-allergenic soft fiber, which provides perfect air permeability, excellent absorption and moisture removal from the body.

Contact us

Our specialists will help you choose the perfect orthopedic pillow for sleeping, adapted to the individual needs of each patient.

ORTHOPEDIC PILLOW – MORFEUSZ CLASSIC 660

- made of elastic foam
- adjustable height
- provides support and correct position of the cervical spine and head during sleep
- degree of hardness: hard
- dimensions: 48 × 30 × 11 cm, with a pad: 48 × 30 × 14 cm



Intended use:

Pillow intended for use in the prophylaxis and prevention of degeneration and pain in the cervical spine caused by an incorrect sleeping position; for those who have problems resulting from the incorrect position of the neck and head during sleep (e.g. back pain, muscle tension, tingling and numbness of the limbs, etc.); who need stable support of the head and cervical spine during sleep, while reducing the pressure on the body.

ORTHOPEDIC PILLOW – MORFEUSZ COMFORT 661

- made of thermoelastic "memory" foam
- adjustable height
- provides support and correct position of the cervical spine and head during sleep
- degree of hardness: medium
- dimensions: 48 × 30 × 11 cm, with a pad: 48 × 30 × 14 cm



Intended use:

Pillow intended for use in the prophylaxis and prevention of degeneration and pain in the cervical spine caused by an incorrect sleeping position; for those who have problems resulting from the incorrect position of the neck and head during sleep (e.g. back pain, muscle tension, tingling and numbness of the limbs, etc.); who need stable support of the head and cervical spine during sleep, while reducing the pressure on the body.

ORTHOPEDIC PILLOW – MORFEUSZ EXCLUSIVE 662

- made of a combination of elastic foam and thermoelastic “memory” foam
- adjustable height
- provides support and correct position of the cervical spine and head during sleep
- degree of hardness: soft
- dimensions: 46 × 30 × 10 cm, with a pad: 46 × 30 × 13 cm

Intended use:

Pillow intended for use in the prophylaxis and prevention of degeneration and pain in the cervical spine caused by an incorrect sleeping position; for those who have problems resulting from the incorrect position of the neck and head during sleep (e.g. back pain, muscle tension, tingling and numbness of the limbs, etc.); who need stable support of the head and cervical spine during sleep, while reducing the pressure on the body.



ORTHOPEDIC PILLOW – MORFEUSZ TRAVEL 663

- made of thermoelastic “memory” foam
- adjustable height
- provides support and correct position of the cervical spine and head during sleep
- perfect for people who often sleep away from home
- degree of hardness: soft
- dimensions: 32 × 25 × 10 cm, with a pad: 32 × 25 × 13 cm

Intended use:

Pillow intended for use in the prophylaxis and prevention of degeneration and pain in the cervical spine caused by an incorrect sleeping position; for those who have problems resulting from the incorrect position of the neck and head during sleep (e.g. back pain, muscle tension, tingling and numbness of the limbs, etc.); who need stable support of the head and cervical spine during sleep, while reducing the pressure on the body.



ORTHOPEDIC PILLOW – MORFEUSZ

SOFT 664

- made of thermoelastic “memory” foam
- adjustable height
- provides support and correct position of the cervical spine and head during sleep
- degree of hardness: soft
- dimensions: 48 × 32 × 11 cm, with a pad: 48 × 32 × 14 cm

Intended use:

Pillow intended for use in the prophylaxis and prevention of degeneration and pain in the cervical spine caused by an incorrect sleeping position; for those who have problems resulting from the incorrect position of the neck and head during sleep (e.g. back pain, muscle tension, tingling and numbness of the limbs, etc.); who need stable support of the head and cervical spine during sleep, while reducing the pressure on the body.



ORTHOPEDIC PILLOW – MORFEUSZ

FOR SLEEP APNEA THERAPY 665

- made of thermoelastic “memory” foam
- adjustable height
- thanks to a specialized profile it prevents the CPAP mask from moving, and thus allows the user to use the device properly
- the special shoulder cut-out allows the user to take a comfortable position, and the raised collar stabilizes the position of the head
- degree of hardness: soft
- dimensions: 55 × 36 × 11 cm, with a pad: 55 × 36 × 14 cm

Intended use:

Orthopedic pillow made of thermoelastic foam, intended for the treatment of sleep apnea, thanks to its special profile it prevents a CPAP mask from moving.







**ANTI-BEDSORE SEAT
CUSHIONS**

Manufacturer of anti-bedsore seat cushions

A sedentary lifestyle is a severe problem of our times. Lack of physical movement is a major cause of many civilization's diseases. Years spent behind a desk can lead to serious problems, including blood circulation and spine issues. Incorrect sitting position is especially dangerous for our health. A very common occurrence is involuntary slouching or sliding down while sitting in an armchair.

How to minimize the negative effect of a sedentary lifestyle?

Anti-bedsore seat cushions help reduce pressure on soft tissues and ensure proper blood flow. They stabilize the pelvis, making it easier to maintain the correct posture.

Which anti-bedsore cushion for sitting to choose – pneumatic or foam?

For people at high risk of developing bedsores, the best solution seems to be a pneumatic cushion. Such cushions are recommended for people who spend most of the day in a sitting position. It is also a good solution for people who have already experienced negative effects of prolonged sitting, but are still able to move independently.

If a patient requires stabilization and minimized pressure on soft tissues, a foam cushion is worth choosing. Such a product adapts better to the shape of the pelvis and is not susceptible to uncontrolled movements like in the case of the pneumatic cushion. It allows to maintain the correct sitting position. Cushions with a gel insert are also a very good solution, as they not only stabilize the pelvis well, but thanks to the gel pad, they redistribute the pressure on the pelvis better, making the user feel more comfortable.

ANTI-BEDSORE SEAT CUSHION BIOFLOTE 2 PLUS

BioFlote™

**Available
colors:**

● black

- high quality pneumatic cushion
- protects users against tissue deformation and bedsores
- made of comfortable, flexible and very durable neoprene
- up to 65% pressure reduction

TECHNICAL PARAMETERS		
MODELS		
reference number	BioFlote™ 2 Plus/40	BioFlote™ 2 Plus/45
dimensions	40 cm × 40 cm × 7 cm	45 cm × 40 cm × 7 cm
max. safe load	210 kg	210 kg



Intended use:

The BioFlote™ 2 Plus anti-bedsores seat cushion protects the user against bedsores and tissue deformation. The cushion adapts to the individual shape of the user. It is intended for wheelchair users, people who stay in a sitting position for a long time, people at risk of developing bedsores or those requiring prophylaxis of pelvic inclination.

ANTI-BEDSORE SEAT CUSHION BIOFLOTE 3 PLUS

BioFlote™

**Available
colors:**

● brown

- reduces the risk of bedsores
- reduces the surface pressure on the pelvis and lumbar spine
- supports the correct sitting position
- universal, convenient and easy to use
- automatically inflates without the use of a pump
- air release valve makes the pillow adapts itself to the shape of the pelvis and adjusts the softness of the seat surface

TECHNICAL PARAMETERS		
MODELS		
reference number	BioFlote™ 3 Plus/43	BioFlote™ 3 Plus/45.5
dimensions	40.5 cm × 43 cm × 7.5 cm	40.5 cm × 45.5 cm × 7.5 cm
max. safe load	130 kg	130 kg



Intended use:

The BioFlote™ 3 Plus anti-bedsores seat cushion is intended for anti-bedsores prevention in people who stay in a sitting position for a long time.

ANTI-BEDSORE RELIEF SEAT CUSHION

HF2008

Available colors:

● black

- chambers of the cushion automatically adjust to the shape of the hips
- reduces pressure on the pelvic area
- disperses the pressure points of the thighs and hips
- provides optimal stability and support
- helps keep blood flowing freely
- improves seating comfort
- protects users from tissue deformation and bedsores
- maximum safe load up to 150 kg



TECHNICAL PARAMETERS

dimensions	42 cm × 42 cm × 5 cm
max. user weight	150 kg

Intended use:

The cushion was designed to fit the shape of the user's body better and reduce pressure. It is intended for patients who are likely to spend more than 15 hours a day in bed and patients at medium to high risk of developing bedsores.

3-LAYER PNEUMATIC SEAT CUSHION

Available colors:

● blue

- easy-to-use pneumatic cushion
- relieves and provides comfort
- unique air channel technology for optimal stability and support
- reduces pressure and distributes weight evenly
- maximum safe load up to 130 kg



TECHNICAL PARAMETERS

dimensions	40 cm (W) × 40 cm (L) × 5.5 cm (H)
weight	198 g
max. user weight	130 kg

Intended use:

The three-layer pneumatic seat cushion is intended for people who spend a lot of time in a sitting position, especially for people at risk of developing bedsores.

ANTI BEDSORE SEAT CUSHION

SUPERŽEL 681

- profiled surface of the cushion ensures optimal distribution of body pressure and facilitates air circulation
- intended for people who spend a lot of time in a sitting position, e.g. at the computer, in a wheelchair, in a car, at work
- perfectly adapts to the pelvis shape
- has a liquid-repellent and vapor-permeable cover with a pocket for small items

TECHNICAL PARAMETERS	
cover	waterproof polyester
dimensions	40 × 42 × 7 cm
	44 × 42 × 7 cm
	48 × 42 × 7 cm



Intended use:

The Super Žel anti-bedsore seat cushion is recommended for anti-bedsore prevention. The unique medical profile and the open structure of the foam allows free air circulation, dissipates heat and moisture on the patient's skin. The gel insert provides pressure redistribution over a wide surface, thus reducing pressure on the body. The bottom layer of the foam dissipates the pressure of the surface on the body and stabilizes the cushion.

ANTI BEDSORE SEAT CUSHION

JUSTYNA 624

- contoured surface stimulates microcirculation
- foam structure reduces the pressure around the coccyx and perineum
- provides excellent air circulation
- eliminates pain and discomfort during prolonged sitting

TECHNICAL PARAMETERS	
cover	100% polyester
dimensions	40 × 42 × 6 cm




Intended use:

The Justyna anti-bedsore seat cushion is recommended in anti-bedsore prevention. It ensures proper air circulation. It is also suitable for people who stay in a sitting position for a long time.

ANTI BEDSORE SEAT CUSHION

SUPERJEŽ 626

- profiled surface of the cushion ensures optimal distribution of body pressure and facilitates air circulation
- skin massaging nubs stimulate blood circulation and prevent bedsore formation
- removable center allows to reduce pressure in the area of the coccyx, anus and reproductive organs
- optionally cotton-polyester or waterproof polyester cover to choose from



Intended use:

This anti-bedsore seat cushion is recommended in anti-bedsore prevention. It ensures proper air circulation and, if necessary, relieves pressure around the intimate body parts. It is also suitable for people who stay in a sitting position for a long time.





ANTI-BEDSORE AND POSITIONING PADS

Manufacturer of anti-bedsore and positioning pads

In hospital care, patients often require area support for anti-bedsore prevention in places particularly exposed to pressure ulcers or relief and stabilization of places where they have already developed. Anti-bedsore and positioning foam pads are an ideal solution. They are easy to use, effective, and their price is more economical than using full anti-bedsore systems, such as anti-bedsore mattresses.

In addition, they support patient positioning, for example, allowing to maintain appropriate abduction of the lower limbs after hip surgery, keep the knees bent or stabilize the patient in a lateral position.

Our offer includes various types of anti-bedsore pads and body positioners such as:

- heel and elbow protectors
- ring positioners
- head positioning pads
- anti-rotation stabilizers
- wedges
- positioning rolls and abduction pillows

HEEL PROTECTOR

ELLEN AD-1

- intended for people at risk of developing pressure ulcers on the heels
- made of polyurethane foam with a characteristic contoured structure
- hook and loop fastening
- protects the foot against abrasions around the heels and ankles
- supports prophylaxis of pressure ulcers
- sold as a set – 2 pcs
- dim. 16 × 19 cm

Intended use:

Foam protector. Provides the reduction of surface pressure on soft tissues around the heels.



ELBOW PROTECTOR

LIMA AD-2

- intended for people at risk of developing pressure ulcers on the elbows
- made of polyurethane foam with a characteristic contoured structure
- hook and loop fastening
- supports prophylaxis of pressure ulcers
- sold as a set – 2 pcs
- dim. 15 × 39.5 cm

Intended use:

Foam protector. Provides the reduction of surface pressure on soft tissues around the elbows.



RING POSITIONER

KARI AD-3

- anti-bedsore ring, in a cover made of soft terry knitted fabric, with hook and loop fastening
- made of polyurethane foam of high density and hardness
- intended for people at risk of developing pressure ulcers on the heels, elbows, head or coccyx
- protects against abrasions, supporting prophylaxis of pressure ulcers
- product sold individually

TECHNICAL PARAMETERS		
dimensions	AD-3/1	Ø 13 cm, height 5 cm
	AD-3/2	Ø 15 cm, height 5 cm
	AD-3/3	Ø 18 cm, height 5 cm
	AD-3/4	Ø 20 cm, height 6 cm
	AD-3/5	Ø 25 cm, height 7 cm
	AD-3/6	Ø 40 cm, height 5 cm
	AD-3/7	Ø 40 cm, height 8 cm

Intended use:

Anti-bedsore foam ring in a cover. It provides adequate stabilization and reduction of surface pressure on soft tissues, around the elbows, heels, head and coccyx.



LEG POSITIONER

ASTI AD-4

- leg positioner is intended for people who require proper and stable positioning of the lower limb, after treatments, operations, paralysis
- made of polyurethane foam with characteristic openings facilitating air circulation
- protects against abrasions, supporting prevention of pressure ulcers, and stabilizes the position of the foot
- dim. 35 × 20 × 10 cm



Intended use:

Anti-bedsore foam positioner. Provides adequate stabilization and reduction of surface pressure on soft tissues around the calves, heels and feet.

HEAD POSITIONER

RUNA AD-5

- head stabilizing positioner
- intended for people at risk of developing pressure ulcers on the head and/or in need of its stabilization
- protects against abrasions, supporting prevention of pressure ulcers, and stabilizes the head in various positions
- made of polyurethane foam, with a characteristic cut-out facilitating air circulation



TECHNICAL PARAMETERS		
dimensions	AD-5/M	18 × 16 × 9 cm
	AD-5/L	23 × 16 × 11 cm

Intended use:

Anti-bedsore foam positioner. It provides adequate stabilization and reduction of surface pressure on soft tissues around the head.

UNIVERSAL ANTI-BEDSORE PAD

MONA AD-6

- universal anti-bedsore pad
- intended for people at risk of developing pressure ulcers, as a base for upper and lower limbs
- protects against abrasions, supporting prophylaxis of pressure ulcers
- made of polyurethane foam with a characteristic contoured structure
- available in a version with and without a cotton-polyester cover
- sold as a set – 2 pcs
- dim. 3 × 20 × 41 cm



Intended use:

Foam anti-bedsore pad, optionally with a cotton-polyester cover. It provides adequate stabilization and reduction of the surface pressure on soft tissues.

REHABILITATION ROLLER

UNNI AD-7

- roller for physical rehabilitation
- works well during various types of rehabilitation and physiotherapeutic treatments, corrective gymnastics, traditional and sports massages
- can be used to stabilize the patient's body, as well as to lift and separate limbs
- in a waterproof washable cover made of artificial leather with a zipper
- additional stitching on the edge, protecting against tearing
- product color sent randomly

TECHNICAL PARAMETERS		
dimensions	AD-7/1	dia. 10 cm, L 40 cm
	AD-7/2	dia. 10 cm, L 50 cm
	AD-7/3	dia. 13 cm, L 40 cm
	AD-7/4	dia. 13 cm, L 50 cm
	AD-7/5	dia. 13 cm, L 60 cm
	AD-7/6	dia. 15 cm, L 50 cm
	AD-7/7	dia. 15 cm, L 60 cm
	AD-7/8	dia. 20 cm, L 60 cm
	AD-7/9	dia. 20 cm, L 100 cm
	AD-7/10	dia. 30 cm, L 100 cm
	AD-7/11	dia. 40 cm, L 100 cm
	AD-7/12	dia. 50 cm, L 100 cm

Intended use:

Foam rehabilitation roller in a washable waterproof cover. Intended for use during various types of rehabilitation treatments, corrective gymnastics or during massage. For use in rehabilitation offices and for home use.



REHABILITATION HALF-ROLLER


AD-8

- half-roller for physical rehabilitation
- works well during various types of rehabilitation and physiotherapeutic treatments, corrective gymnastics, traditional and sports massages
- can be used to stabilize the patient's body, as well as to lift and separate limbs
- in a waterproof washable cover made of artificial leather with a zipper
- additional stitching on the edge, protecting against tearing
- product color sent randomly

TECHNICAL PARAMETERS		
dimensions	AD-8/1	dia. 12/6 cm, L 40 cm
	AD-8/2	dia. 15/7.5 cm, L 50 cm
	AD-8/3	dia. 15/7.5 cm, L 60 cm
	AD-8/4	dia. 20/10 cm, L 60 cm
	AD-8/5	dia. 40/20 cm, L 60 cm

Intended use:

Foam rehabilitation half-roller in a washable waterproof cover. Intended for use during various types of rehabilitation treatments, corrective gymnastics or during massage. For use in rehabilitation offices and for home use.



TRIANGULAR POSITIONING PAD

PP-8

- triangular anti-bedsore pad
- intended for people who require proper body position stabilization (e.g. bent knees, lateral body position)
- stabilizes the relevant body parts, supporting anti-bedsore prophylaxis
- made of polyurethane foam
- available with or without a cotton-polyester cover
- dim. 19 × 19 × 43 cm




Intended use:

Triangular foam anti-bedsore pad, intended to stabilize the body position, e.g. bent knees, lateral body position or abduction of the upper or lower limb.

30°/60° WEDGE

PP-9

- 30°/60° wedge
- intended for people who require a higher torso position in relation to the rest of the body or who stay in bed in a sitting position
- provides comfortable and safe back support
- used for massage, rehabilitation or corrective gymnastics
- made of polyurethane foam
- available with or without a cotton-polyester cover
- dim. 22.5 × 60 × L 58 cm



Intended use:

Anti-bedsore foam 60° wedge intended for people who require a higher torso position in relation to the rest of the body or stabilization in a specific position, e.g. bent knees or lateral position.

45° WEDGE

PP-10

- 45° wedge
- intended for people who require a higher torso position in relation to the rest of the body or stabilization (e.g. bent knees, lateral body position)
- provides comfortable stabilization and safe support for the back, allows for higher body position
- used for massage, rehabilitation or corrective gymnastics
- made of polyurethane foam
- available with or without a cotton-polyester cover
- dim. 44.5 × 44.5 × 60 cm



Intended use:

Anti-bedsore foam 45° wedge intended for people who require a higher torso position in relation to the rest of the body or stabilization in a specific position, e.g. bent knees or lateral position.

LONG LEG ABDUCTION PILLOW

PP-11

- intended for people who require immobilization and stabilization of legs in relation to each other and keeping them in a parted position, abduction of lower limbs after surgery (e.g. endoprosthesis of the hip joint)
- protects against abduction, adduction and rotation movements in the hip joints, in the postoperative period after lower limb injuries and supports the legs along their entire length in the right position, without causing pressure
- made of polyurethane foam with soft straps with hook and loop fastening

TECHNICAL PARAMETERS		
dimensions	PP-11/S	W 31 × height 46 cm
	PP-11/M	W 36 × height 56 cm
	PP-11/L	W 41 × height 64 cm



Intended use:

Long foam leg abduction pillow with elastic straps with hook and loop fastening, intended for people who require immobilization and stabilization of legs in relation to each other and keeping them in a parted position, abduction of lower limbs after surgery (e.g. hip joint endoprosthesis).

SHORT LEG ABDUCTION PILLOW

PP-12

- short leg abduction pillow is intended for people who require immobilization and stabilization of legs in relation to each other and keeping them in a parted position, abduction of lower limbs after surgery (e.g. hip joint endoprosthesis)
- supports legs in the right position without causing pressure
- made of polyurethane foam, in a cover made of terry knitted fabric, with hook and loop fastener

TECHNICAL PARAMETERS		
dimensions	PP-11/M	17.5 × 14 × 24 cm
	PP-11/L	20 × 16.5 × 26 cm



Intended use:

Short foam leg abduction pillow with hook and loop fastener, intended for people who require immobilization and stabilization of legs in relation to each other and keeping them in a parted position, abduction of lower limbs after surgery (e.g. hip joint endoprosthesis), for pregnant women or in the postpartum period.





WALKING AIDS

Manufacturer of walking frames, handcarts, walkers, crutches and canes

Often, during the rehabilitation period after an injury or in case of instability while walking or standing, the patient looks for additional support that will assist him when moving or standing up.

In our wide offer our customers will find rehabilitation walking frames, walking frames for disabled people, walking frames and walkers for seniors as well as orthopedic crutches supporting walking, canes, tripods, tetrapods and pedal exercisers at an affordable price.

Rehabilitation walking frames for adults, disabled people and seniors

It is an excellent solution to support movement mainly around the house. Walking frames can be divided into 3 main categories:

- walking frames that need to be lifted up for movement – the best support for a standing position, however, it requires being picked up from the ground and moved forward as the user walks, which at the same time causes restrictions in movement.
- walking frames that move with you – perfect for seniors or adults learning to walk, e.g. after strokes. It easily moves along with the user, providing support all the time.
- walking frames with wheels – perfect for seniors or people with disabilities who are able to move quite confidently. These rehabilitation walking frames are equipped with two front wheels to facilitate patient movement.

Forearm walking frames

They are used to bring support when walking and to make it easier for the user to move around without the help of another person. They provide support all over the patient's arms, which significantly relieves the lower limbs.

Handcarts for walking

This solution combines the features of walking frames and walkers. Handcarts have two legs, two large wheels and a seat, which enable moving around the house as well as outdoors for short distances.

Orthopedic crutches for walking

Rehabilitation crutches are primarily used as support for people who need to relieve the lower limbs, e.g. after fractures or sprains or additional stabilization while moving. Their advantage is a much smaller size than walking frames or walkers, however, they are not self-standing products, which means less support in comparison to walkers or walking frames. Orthopedic crutches can be divided into two types: elbow crutches, where the weight of the body rests mostly on the hands, and arm crutches, where some of the weight is transferred to the armpits.

Tripod and Tetrapod

A special variety of canes are tripods and tetrapods offering more points of support, which allow the patient to be better supported than with a standard cane.

Walking canes

It is an ideal support for seniors and people in the final period of rehabilitation. It helps maintain balance and partially relieves the lower limb. Walking canes are exceptionally light and easy to carry. In our offer you can also find foldable canes.

WALKING FRAME

RF-100

- lightweight aluminum walking frame for people with standing problems
- helps with daily activities, e.g. getting out of bed, walking, using the toilet, taking a stroll
- folded by pulling up a wooden ball
- ergonomic plastic handles
- equipped with two front wheels
- adjustable height

TECHNICAL PARAMETERS	
max. safe load	127 kg
height adjustment	78.7-86.3 cm
maximum width	64 cm
weight	1.94 kg
depth	51 cm



Available colors:

- gray
- blue
- brown

Intended use:

The walking frame is intended for elderly and disabled people with diseases of the locomotor system or with the problem of verticalization. It can be used inside. It helps with daily activities, e.g. getting out of bed, walking, using the toilet.

WALKING FRAME

RF-131

- lightweight aluminum walking frame for people with standing problems
- helps with daily activities, e.g. getting out of bed, walking, using the toilet, taking a stroll
- combines 3 functions – sliding (moves with the user), fixed (has a rigid frame), transport (option with and without wheels)
- the sliding and fixed function can be conveniently changed using the locking button
- soft foam handles
- adjustable height

TECHNICAL PARAMETERS	
reference number	RF-131
height (without / with wheels)	77-95 / 79-97 cm
total width (without wheels / with wheels)	57 / 60 cm
total length (without wheels / with wheels)	46 / 49 cm
front wheel diameter (without wheels / with wheels)	N/A / 5 inches
max. safe load (without wheels / with wheels)	100 kg / 100 kg
weight (without wheels / with wheels)	2.4 kg / 2.75 kg



Available colors:

- gray
- blue
- brown

Intended use:

The walking frame is intended for elderly and disabled people with diseases of the locomotor system or with the problem of verticalization. It can be used inside. It helps with daily activities, e.g. getting out of bed, walking, using the toilet.

WALKING FRAME WITH BENCH AND TOILET SET – HARRY 133

RF-133

- steel-aluminum walking frame for people with standing problems
- helps with daily activities
- comes with a bench and a toilet set included, so it can also be used as a toilet seat and a shower chair at the same time
- handles placed at two heights, providing additional help when getting up and sitting down
- adjustable height

TECHNICAL PARAMETERS	
height adjustment	81-99 cm
frame width	52 cm
toilet set width	40 cm
bench width	41 cm
seat height	46 cm
max. safe load	100 kg
weight	5.9 kg



Available colors:

- gray

Intended use:

The walking frame is intended for elderly and disabled people with diseases of the locomotor system or with the problem of verticalization. It can be used inside. It helps with everyday activities, e.g. getting out of bed, walking, using the toilet.

WALKING FRAME – LUCA 803/1

RF-803/1

- steel-aluminum walking frame for people with standing problems
- helps with daily activities
- handles placed at two heights, providing additional help when getting up and sitting down
- adjustable height
- max. user weight 100 kg

TECHNICAL PARAMETERS	
height adjustment	81-99 cm
frame width	52 cm
max. safe load	100 kg
weight without wheels	2.5 kg
weight with wheels	3.5 kg

Intended use:

The walking frame is intended for elderly and disabled people with diseases of the locomotor system or with the problem of verticalization. It can be used inside. It helps with everyday activities, e.g. getting out of bed, walking, using the toilet.

Available colors:

● gray



FOREARM WALKING FRAME – ATLANTA 141

RF-141

- steel forearm walking frame for people with standing problems
- used also to aid walking with an assistant and to increase the user's ability to move around without the presence of another person
- brakes on the rear wheels
- soft foam handles
- drip holder

TECHNICAL PARAMETERS	
max. safe load	100 kg
width	64 cm
width when folded	21 cm
length	73 cm
height adjustment	113-135.5 cm
weight	10.8 kg
heel diameter	10 cm

Intended use:

The forearm walking frame is intended for the elderly and disabled people with diseases of the locomotor system or as an aid during the verticalization process. It is used to aid walking and to increase the user's ability to move without the presence of another person. It can be used indoors.

Available colors:

● gray



2-WHEELED HANDCART WITH SEAT

RF-651

- steel handcart for people with standing problems
- equipped with a seat so that the user can rest on
- ergonomic plastic handles
- two solid wheels at the front

TECHNICAL PARAMETERS	
max. safe load	113 kg
handle height adjustment	78-98 cm
total width	55 cm
depth	52 cm
weight	5 kg

Intended use:

The handcart is intended for elderly and disabled people with diseases of the locomotor system or with the problem of verticalization. It helps people who need support while walking, and in the event of fatigue, it also gives the opportunity to rest. It can be used inside. It helps with everyday activities, e.g. getting out of bed, walking, using the toilet.

Available colors:

● gray



4-WHEELED WALKER

RF-610

- steel 4-wheeled walker
- solid and stable frame
- with a soft seat
- with a convenient bag under the seat
- detachable foam backrest
- with a pocket for small items under the seat

TECHNICAL PARAMETERS	
max. safe load	136 kg
handle height adjustment	79-92.5 cm
width when unfolded	61.5 cm
width when folded	27.5 cm
total length	64 cm
wheel diameter	20 cm
weight	7.8 kg
width between handles	450 mm
width between rear wheels	530 mm
seat height	515 mm
seat depth	340 mm
seat width	300 mm
width when folded	615 mm
length when folded	175 mm
load capacity of the bag	5 kg

Intended use:

The walker is intended to help people with limited mobility who require support when getting up, standing and walking. It enables independent performance of most daily motor activities. It can be used indoors and outdoors.

Available colors:

- gray
- red
- blue
- black



4-WHEELED WALKER – DENVER PREMIUM 615

RF-615

- aluminum 4-wheeled walker
- maneuverable frame
- with a soft seat
- with a portable basket under the seat, which can also be mounted on the front of the frame
- detachable backrest finished with foam

TECHNICAL PARAMETERS	
max. safe load	136 kg
handle height adjustment	82-93 cm
width when unfolded	61 cm
width when folded	23 cm
total length	69 cm
wheel diameter	20 cm
weight	7.5 kg
folded length	99 cm
width between handles	476 mm
width between rear wheels	460 mm
seat height	575 mm
load capacity of the basket	5 kg
seat width	355 mm
seat depth	340 mm

Intended use:

The walker is intended to help people with limited mobility who require support when getting up, standing and walking. It enables independent performance of most daily motor activities. It can be used indoors and outdoors.

Available colors:

- gray
- red



4-WHEELED WALKER – OSLO PREMIUM 616 RF-616

- aluminum 4-wheeled walker
- maneuverable frame
- with a soft seat
- with a portable basket under the seat, which can also be mounted on the front of the frame
- detachable backrest finished with foam
- backrest can be installed in both directions
- flip-up footrest, which is used to transport the user

TECHNICAL PARAMETERS	
max. safe load	136 kg
handle height adjustment	81-91 cm
width when unfolded	65 cm
width after unfolding the footrest	88 cm
width when folded	23 cm
total length	66 cm
wheel diameter	20 cm
weight	8.2 kg
folded length	99 cm
width between handles	476 mm
width between rear wheels	460 mm
seat height	575 mm
load capacity of the basket	5 kg
seat width	355 mm
seat depth	340 mm

Intended use:

The walker is intended to help people with limited mobility who require support when getting up, standing and walking. It enables independent performance of most daily motor activities. It can be used indoors and outdoors.

Available colors:

- gray
- red



4-WHEELED WALKER – LIZBONA PREMIUM 620 RF-620

- high-quality 4-wheeled walker
- modern and sleek design
- lightweight and maneuverable
- easy to move around
- foldable aluminum frame
- foldable to small size
- comfortable seat
- with brakes
- with a large shopping bag

TECHNICAL PARAMETERS	
max. safe load	136 kg
handle height adjustment	81-93 cm
width when unfolded	67 cm
total length	70 cm
wheel diameter	20 cm
weight	6.9 kg
seat height	54.5 cm
seat width	45 cm
width between handles	51.5 cm
width between rear wheels	58 cm
seat height	54.5 cm
seat width	45 cm
seat depth	25 cm
width when folded	23 cm
load capacity of the bag	5 kg

Intended use:

The walker is intended to help people with limited mobility who require support when getting up, standing and walking. It enables independent performance of most daily motor activities. It can be used indoors and outdoors.

Available colors:

- gray
- red



4-WHEELED WALKER – PALERMO PREMIUM 624

RF-624

Available colors:

- gray
- blue

- economical 4-wheeled walker
- modern and sleek design
- foldable aluminum frame
- easy to transport
- foldable to a very small size
- comfortable seat
- with brakes
- with a handy bag

TECHNICAL PARAMETERS	
max. safe load	136 kg
handle height adjustment	80-96 cm
width when unfolded	64.5 cm
wheel diameter	20 cm
seat height	53 cm
seat depth	25 cm
seat width	46 cm
weight	8.2 kg
total length	66 cm
total width	64.5 cm
length when folded	64 cm
width when folded	27 cm
height when folded	34 cm

Intended use:

The walker is intended to help people with limited mobility who require support when getting up, standing and walking. It enables independent performance of most daily motor activities. It can be used indoors and outdoors.



4-WHEELED WALKER – BOGOTA PREMIUM 625

RF-625

Available colors:

- gray

- 4-wheeled walker
- modern and sleek design
- foldable aluminum frame
- reversible backrest
- with shopping bag and crutch holder
- swing away and detachable footrests
- transport function

TECHNICAL PARAMETERS	
max. safe load	136 kg
handle height adjustment	85-98 cm
overall width	640 mm
weight	9.3 kg
total length	74 cm
total width	62 cm
width when folded	30 cm
height when folded	85 cm
wheel diameter	20 cm
seat height	53 cm
seat depth	25 cm
seat width	46 cm

Intended use:

The walker is intended to help people who need support while walking, and in the event of fatigue, it also provides the opportunity to rest. It can be used indoors and outdoors. The walker can also be used to transport the user by an assistant.



3-WHEELED WALKER – BARCELONA 611

RF-611

- stable steel foldable frame
- with a bag, basket and brakes
- PVC + PP wheels (polyvinyl chloride + polypropylene)

TECHNICAL PARAMETERS	
max. safe load	110 kg
handle height adjustment	82-97 cm
width when unfolded	60 cm
width when folded	25 cm
total length	62 cm
wheel diameter	20 cm
weight	7.4 kg

Intended use:

The walker is intended to help people with limited mobility who require support when getting up, standing and walking. It enables independent performance of most daily motor activities. It can be used indoors and outdoors.



Available colors:

- gray
- blue

3-WHEELED WALKER – HELSINKI PREMIUM 621

RF-621

- lightweight aluminum foldable frame
- with bag, basket and brakes
- PU + PP wheels (polyurethane + polypropylene)

TECHNICAL PARAMETERS	
max. safe load	110 kg
handle height adjustment	82-97 cm
width when unfolded	56 cm
width when folded	25 cm
total length	62 cm
wheel diameter	20 cm
weight	5.8 kg

Intended use:

The walker is intended to help people with limited mobility who require support when getting up, standing and walking. It enables independent performance of most daily motor activities. It can be used indoors and outdoors.



Available colors:

- blue

3-WHEELED WALKER – GENEWA PREMIUM 622

RF-622

- modern and sleek design
- lightweight aluminum foldable frame
- with a bag and brakes
- handle height adjustment with locking buttons
- TPR wheels (thermoplastic rubber)

TECHNICAL PARAMETERS	
max. safe load	115 kg
handle height adjustment	82.5-92.5 cm
total width	69.5 cm
total length	66 cm
weight	5 kg
wheel diameter	20 cm

Intended use:

The walker is intended to help people with limited mobility who require support when getting up, standing and walking. It enables independent performance of most daily motor activities. It can be used indoors and outdoors.



Available colors:

- gray
- brown

ELBOW CRUTCH

RF-200

- hinged cuff
- lower height adjustment

TECHNICAL PARAMETERS	
max. safe load	180 kg
height adjustment	66-91.5 cm
weight	0.53 kg
height from the handle to the cuff	22 cm

Intended use:

Elbow crutches are intended for movement, belaying while walking, for supporting oneself, maintaining balance and transferring body weight in cases of traumatic and orthopedic diseases of the lower limbs when full load on the limb is prohibited. Intended for patients after injuries, accidents, fractures, sprains, dislocations, for seniors.



Available colors:

- gray

ELBOW CRUTCH

RF-210

- hinged cuff
- lower and upper height adjustment

TECHNICAL PARAMETERS	
max. safe load	180 kg
height adjustment above the handle	26.5-34 cm
height adjustment below the handle	63.5-91 cm
weight	0.55 kg

Intended use:

Elbow crutches are intended for movement, belaying while walking, for supporting oneself, maintaining balance and transferring body weight in cases of traumatic and orthopedic diseases of the lower limbs when full load on the limb is prohibited. Intended for patients after injuries, accidents, fractures, sprains, dislocations, for seniors.



Available colors:

- gray

ELBOW CRUTCH

RF-220

- open cuff
- lower height adjustment
- wide range of colors

TECHNICAL PARAMETERS	
max. safe load	113 kg
height adjustment	73.6-96.5 cm
weight	0.62 kg

Intended use:

Elbow crutches are intended for movement, belaying while walking, for supporting oneself, maintaining balance and transferring body weight in cases of traumatic and orthopedic diseases of the lower limbs when full load on the limb is prohibited. Intended for patients after injuries, accidents, fractures, sprains, dislocations, for seniors.



Available colors:

- gray
- blue
- red
- green
- pink
- black

UNDERARM CRUTCH

RF-230

- underarm (armpit) crutch with high adjustment
- sizes M, L

TECHNICAL PARAMETERS		
reference number	RF-230/M	RF-230/L
max. safe load	113 kg	113 kg
weight	1.5 kg	1.7 kg
height adjustment	114-134 cm	134-155 cm

Intended use:

Underarm crutches are intended for movement, belaying while walking, for supporting oneself, maintaining balance and transferring body weight in cases of traumatic and orthopedic diseases of the lower limbs when full load on the limb is prohibited. Intended for patients after injuries within the musculoskeletal system, after surgeries, fractures, with paresis of the lower limbs.



Available colors:

- gray

TRIPOD

RF-300

- adjustable height
- with a wide stable leg spacing
- plastic ergonomic handle

TECHNICAL PARAMETERS	
max. safe load	127 kg
height adjustment	76-01 cm
weight	1.15 kg
distance between leg tips	30 cm; 36 cm

Intended use:

The tripod is intended for movement, belaying while walking, for supporting, maintaining balance and transferring body weight in cases of traumatic and orthopedic diseases of the lower limbs. The tripod has three support points, which provide the user with greater stability while walking. It is intended for patients with balance problems, after injuries, during the treatment of degenerative changes, after surgeries, when it is advisable to relieve one of the limbs.



Available colors:

● gray

TETRAPOD

RF-400

- adjustable height
- with a leg spacing that fits the stairs
- plastic ergonomic handle

TECHNICAL PARAMETERS	
max. safe load	127 kg
height adjustment	77-103 cm
weight	1.15 kg
leg height	14.5 cm
distance between leg tips	21 cm

Intended use:

The tetrapod is intended for movement, belaying while walking, for supporting, maintaining balance and transferring body weight in cases of traumatic and orthopedic diseases of the lower limbs. The tetrapod has four support points, which provide the user with greater stability while walking. It is intended for patients with balance problems, after injuries, during the treatment of degenerative changes, after surgeries, when it is advisable to relieve one of the limbs.



Available colors:

● gray

ALUMINUM CANE

RF-510

- with double height adjustment

TECHNICAL PARAMETERS	
max. safe load	113 kg
weight	0.38 kg
height adjustment	78.7-101.6 cm

Intended use:

The cane is intended to help in moving, belay while walking, maintain balance and transfer body weight in cases of traumatic and orthopedic diseases of the lower limbs. It is intended for people with slight disabilities of the locomotor system who need help in moving around.



Available colors:

● black
● blue
● red
● brown
● gray
● green

FOLDABLE ALUMINUM CANE

RF-520

- with double height adjustment

TECHNICAL PARAMETERS	
max. safe load	113 kg
weight	0.35 kg
height adjustment	78.7-89 cm

Intended use:

The cane is intended to help in moving, belay while walking, maintain balance and transfer body weight in cases of traumatic and orthopedic diseases of the lower limbs. It is intended for people with slight disabilities of the locomotor system who need help in moving around.



Available colors:

● black
● blue
● red
● brown
● gray
● green

PEDAL EXERCISER WITH DIGITAL DISPLAY 901

RF-901

- for safe muscle strengthening, muscle conditioning exercises, stimulation of the circulatory system
- intended to exercise upper and lower limbs
- with counter
- counter functions: number of revolutions performed, training duration, estimated number of burned calories, number of revolutions per minute
- with anti-slip tips
- equipped with a resistance mechanism
- foldable frame

TECHNICAL PARAMETERS	
weight	2 kg
width	56.5 cm
height	39.5 cm
length	40 cm
max. user weight	100 kg

Intended use:

The pedal exerciser is intended for active (resistance) exercises of the lower and upper limbs in traumatic, orthopedic and rheumatic diseases. It is used in rehabilitation wards, in nursing homes and for exercises at home. The pedal exerciser is intended for safe muscle strengthening, muscle conditioning exercises, and stimulation of the circulatory system.

Available colors:

gray



PEDAL EXERCISER 902

RF-902

- for safe muscle strengthening, muscle conditioning exercises, stimulation of the circulatory system
- intended to exercise upper and lower limbs
- with anti-slip tips
- equipped with a resistance mechanism

TECHNICAL PARAMETERS	
weight	2 kg
width	56.5 cm
height	39.5 cm
length	40 cm
max. user weight	100 kg

Intended use:

The pedal exerciser is intended for active (resistance) exercises of the lower and upper limbs in traumatic, orthopedic and rheumatic diseases. It is used in rehabilitation wards, in nursing homes and for exercises at home. The pedal exerciser is intended for safe muscle strengthening, muscle conditioning exercises, and stimulation of the circulatory system.

Available colors:

gray



FOLDABLE PEDAL EXERCISER 903

RF-903

- for safe muscle strengthening, muscle conditioning exercises, stimulation of the circulatory system
- intended to exercise upper and lower limbs
- foldable frame
- equipped with a resistance mechanism
- with anti-slip tips
- white sleek frame

TECHNICAL PARAMETERS	
weight	2.1 kg
width	42 cm
height	30 cm
length	23 cm
max. user weight	100 kg

Intended use:

The pedal exerciser is intended for active (resistance) exercises of the lower and upper limbs in traumatic, orthopedic and rheumatic diseases. It is used in rehabilitation wards, in nursing homes and for exercises at home. The pedal exerciser is intended for safe muscle strengthening, muscle conditioning exercises, and stimulation of the circulatory system.

Available colors:

white







WHEELCHAIRS

Manufacturer of wheelchairs

We are convinced that every customer will find a wheelchair in our assortment that will meet their expectations both in terms of functionality and the expected price.

Each of our wheelchairs has been designed to provide the highest possible comfort and functionality to all people suffering from mobility impairments. We took care of a wide range of parameters, thanks to which every user will be able to choose a product adapted to individual needs.

Our wheelchairs are equipped with additional functionalities that affect the comfort of use:

- brake for the assistant
- height adjustable armrests
- backrest tilt angle adjustment
- swing away backrest
- anti-tip wheels
- wheelchair headrest
- lift-up footrests

Wheelchair rental

We meet the needs of patients who require short-term mobile assistance or do not want to pay for the purchase of a wheelchair.

We offer the possibility of using our rehabilitation equipment rental.

The minimum lease period is one month with the possibility of multiple renewals.

Find out more!

Visit our website www.rehafund.pl

STANDARD STEEL WHEELCHAIR

BME4611

Available colors:

● black

- solid or pneumatic wheels
- foldable frame and single cross brace
- plastic handrims
- long armrests
- pocket for small items

TECHNICAL PARAMETERS		
seat width [mm]		450
total width [mm]		620
total length [mm]		1020
total height [mm]		860
seat depth [mm]		420
backrest height [mm]		420
armrest height from the seat [mm]		290
seat height from the ground [mm]		460
max. safe load [kg]		110 kg
basic wheelchair weight [kg]		17 kg
static stability [°]	downhill	9°
	uphill	9°
	transverse	9°
wheel diameter	front	200
	rear	600
use and storage temperature		from -15°C to 40°C



Intended use:

The manual wheelchair is used as a means of transport for the disabled, patients and the elderly.

STEEL WHEELCHAIR

CRUISER 1 BASIC, RF-1-B

Available colors:

● gray

- practical, functional and easy to use steel wheelchair
- foldable frame and single cross brace
- nylon upholstery
- swing away side panels
- swing away and detachable footrests with adjustable high
- steel handrims
- pocket for small items
- solid or pneumatic rear wheels

TECHNICAL PARAMETERS			
seat width [mm]	420	450	480
total width [mm]	605	630	655
total length [mm]	1030	1030	1030
total height [mm]	910	910	910
seat depth [mm]	400	400	400
backrest height [mm]	410	410	410
armrest height from the seat [mm]	250	250	250
max. safe load [kg]	120	120	120
weight in packaging / own weight [kg]	18.1	18.3	18.5
static stability – border value [°]	forward (wheels unlocked)	21.1°	23.3°
	backwards (wheels unlocked)	26.3°	28.2°
	sideways (wheels unlocked)	21.3°	22.7°
wheel diameter [mm]	front	200	200
	rear	610	610
use and storage temperature	-5°C to 55°C	-5°C to 55°C	-5°C to 55°C



Intended use:

The Cruiser 1 Basic wheelchair is intended for people with a permanent disability that limits unassisted walking.

Possibility of choosing additional equipment:

RF-1-B with lift-up footrests

STEEL WHEELCHAIR

CRUISER 1, RF-1

Available colors:

● gray

● black



- steel wheelchair with additional functionalities for the patient
- double cross brace for greater stability
- quick-release rear wheels for easy transportation
- aluminum handrims for easy maneuvering
- foldable frame
- nylon upholstery
- swing away and detachable footrests with adjustable high
- pocket for small items
- solid or pneumatic rear wheels
- foam cushion included

Intended use:

The Cruiser 1 wheelchair is intended for people with a permanent disability that limits unassisted walking.

Possibility of choosing additional equipment:

RF-1-B with flip-up footrests

TECHNICAL PARAMETERS

seat width [mm]		420	450	480	510
total width [mm]		605	630	655	680
total length [mm]		1030	1030	1030	1030
total height [mm]		910	910	910	910
seat depth [mm]		400	400	400	400
backrest height [mm]		410	410	410	410
armrest height from the seat [mm]		250	250	250	250
max. safe load [kg]		120	120	120	135
weight in packaging / own weight [kg]		18.1	18.3	18.5	20
static stability – border value [°]	forward (wheels unlocked)	21.1°	23.3°	19°	19°
	backwards (wheels unlocked)	26.3°	28.2°	26.5°	26.5°
	sideways (wheels unlocked)	21.3°	22.7°	21°	21°
wheel diameter [mm]	front	200	200	200	200
	rear	600	600	600	600
use and storage temperature		-5°C to 55°C	-5°C to 55°C	-5°C to 55°C	-5°C to 55°C

STEEL WHEELCHAIR

CRUISER 2, RF-2

Available
colors:

● black



- with high-quality technical solutions
- single cross brace with a vibration-damping mounting system
- quick-release pins for easy transportation
- swing away and detachable side panels
- weight of 16.5 kg with all accessories, which is an excellent result for a steel wheelchair
- additionally reinforced armrests with a metal insert, which makes them durable
- gravity center adjustable on the front and rear wheels
- height-adjustable footrest
- aluminum handrims
- ergonomic handles for the assistant
- footrests: detachable and can be turned inward and outward
- foam cushion included

Intended use:

The Cruiser 2 wheelchair is intended for people with a permanent disability that limits unassisted walking.

Possibility of choosing additional equipment:

RF-2 with brakes

TECHNICAL PARAMETERS				
seat width [mm]	420	450	480	
total width [mm]	front – 63 cm / rear – 65 cm	front – 66.5 cm / rear – 69 cm	front – 69 cm / rear – 71 cm	
total length [mm]	1110	1110	1110	
total height [mm]	950	950	950	
seat depth [mm]	410	410	410	
backrest height [mm]	450	450	450	
armrest height from the seat [mm]	220	220	220	
max. safe load [kg]	120	120	120	
weight in packaging / own weight [kg]	16.5	16.5	16.5	
wheel diameter [mm]	front	200	200	200
	rear	600	600	600
use and storage temperature	–5°C to 40°C		–5°C to 40°C	

STANDARD ALUMINUM WHEELCHAIR

DY01908LAJ

Available
colors:

● black



- solid or pneumatic wheels
- foldable frame and double cross brace
- metal handrim
- shortened sides that allow for convenient access to the table
- swing away side panels
- pocket for small items

Intended use:

Manual wheelchairs are used as a means of transport for the disabled, the sick and the elderly and infirm people.

TECHNICAL PARAMETERS

total length	97 cm
total width	60, 63, 68 cm
total height	94 cm
seat width	42, 45, 48 cm
front wheel	8 inch
rear wheel	24 inch
seat height	49 cm
seat depth	40 cm
backrest height	41 cm
width when folded	33 cm
net weight	13.7 kg
max. user weight	100 kg

ALUMINUM WHEELCHAIR

CRUISER ACTIVE, RF-3

Available
colors:

● black



- functional, durable and lightweight aluminum wheelchair
- brakes for the assistant
- foldable backrest that helps in transporting the wheelchair
- double cross brace for greater stability
- height-adjustable armrests
- swing away side panels
- detachable and height-adjustable footrests
- quick-release rear wheels
- aluminum handrims
- foam cushion included
- anti-tip wheels
- hip belt

Intended use:**RF-3:**

The Cruiser Active wheelchair is intended for people with permanent disabilities limiting independent walking.

RF-3-Z:

The Cruiser Active Z wheelchair is intended for people with locomotor dysfunctions that prevent walking. Equipping the wheelchair with a headrest provides support for patients requiring head and back stabilization.

Possibility of choosing additional equipment:

RF-3 with a headrest

TECHNICAL PARAMETERS					
seat width [mm]		420	450	480	510
total width [mm]		625	650	675	705
total length [mm]		1030	1030	1030	1030
total height [mm]		950	950	950	950
seat depth [mm]		415	415	415	415
backrest height [mm]		450	450	450	450
armrest height adjustment [mm]		240–310	240–310	240–310	240–310
max. safe load [kg]		135	135	135	135
basic wheelchair weight [kg]		13.3	13.8	14.3	14.8
static stability [°]	downhill	18.5°	19.3°	17.1°	18.1°
	uphill	27.1°	28.6°	26.8°	27.3°
	transverse	18.4°	19.6°	18.2°	18.5°
wheel diameter [mm]	front	200	200	200	200
	rear	610	610	610	610
use and storage temperature		from –5°C to 55°C	from –5°C to 55°C	from –5°C to 55°C	from –5°C to 55°C

ALUMINUM WHEELCHAIR – CRUISER SMART

RF-7

- ingeniously designed wheelchair with an interesting specification
- modern and sleek design
- brakes for the assistant
- swing away side panels
- comfortable armrests with space for hands
- solid or pneumatic wheels
- lightweight and maneuverable

TECHNICAL PARAMETERS		
seat width [mm]	450	
total width [mm]	630	
total length [mm]	1100	
total height [mm]	930	
seat depth [mm]	400	
backrest height [mm]	455	
armrest height from the seat [mm]	270	
seat height from the ground [mm]	535	
max. safe load [kg]	140	
basic wheelchair weight [kg]	14	
static stability [°]	downhill	≥5°
	uphill	≥5°
	transverse	≥5°
wheel diameter	front	200
	rear	600
use and storage temperature	5°C ~ 40°C	

Available colors:

● gray

Upholstery colors:

● blue

● red



Intended use:

The Cruiser Smart wheelchair is intended for people with permanent disabilities limiting independent walking.

ALUMINUM WHEELCHAIR – CRUISER AREO

RF-6

- modern and sleek design
- lightweight, maneuverable and very durable
- double vibration-damping cross brace
- quick-release pin for easy transport in a boot
- weight of the wheelchair with accessories from 14.1 kg
- gravity center adjustable on the front and rear wheels
- height-adjustable footrest
- aluminum handrims
- solid and pneumatic rear wheels
- brakes and ergonomic handles for the assistant
- foam cushion included
- anti-tip wheels

TECHNICAL PARAMETERS			
seat width [mm]	420	450	
total width [mm]	590	640	
total length [mm]	1070	1070	
total height [mm]	890	890	
seat depth [mm]	440	440	
backrest height [mm]	430	430	
armrest height from the seat [mm]	200	200	
seat height from the ground [mm]	520	520	
max. safe load [kg]	130	130	
basic wheelchair weight [kg]	14.1	14.6	
static stability [°]	downhill	≥10°	≥10°
	uphill	≥10°	≥10°
	transverse	≥10°	≥10°
wheel diameter [mm]	front	200	200
	rear	600	600
use and storage temperature	-10°C~40°C		-10°C~40°C

Available colors:

● black



Intended use:

The Cruiser Areo wheelchair is intended for people with permanent disabilities limiting independent walking and maintained functionality in the upper limbs, enabling independent movement in a wheelchair.

ALUMINUM WHEELCHAIR

CRUISER ACTIVE 2, RF-4

Available colors:

- gray
- black



- with high-quality technical solutions
- double cross brace construction with a vibration-damping mounting system
- quick-release pins for easy transport of the wheelchair
- swing away and detachable side panels
- weight of 16 kg with all accessories
- additionally reinforced armrests with a metal insert, which makes them durable
- height-adjustable armrests
- gravity center adjustable on the front and rear wheels
- height-adjustable footrests
- aluminum handrims
- pneumatic and PU rear wheels
- brakes and ergonomic handles for the assistant
- footrests: detachable and can turned outwards and inwards
- foam cushion included
- foldable backrest
- anti-tip wheels
- hip belt

Intended use:

RF-4:

The RF-4 Cruiser Active 2 wheelchair is intended for people with permanent disabilities limiting independent walking and maintained functionality in the upper limbs, enabling independent movement in a wheelchair.

RF-4-Z:

The RF-4 Cruiser Active 2 Z wheelchair is intended for people with locomotor dysfunctions that prevent walking. Equipping the wheelchair with a headrest provides support for patients requiring head and back stabilization.

Possibility of choosing additional equipment:

RF-4 with a headrest

TECHNICAL PARAMETERS			
seat width [mm]	420	450	480
total width [mm]	630	660	690
total length [mm]	1000	1000	1000
total height [mm]	950	950	950
seat depth [mm]	410	410	410
backrest height [mm]	420	420	420
armrest height from the seat [mm]	190–270	190–270	190–270
armrest height adjustment [mm]	80	80	80
max. safe load [kg]	120	120	120
basic wheelchair weight [kg]	16	16	16
ability to climb an obstacle [mm]	60	60	60
wheel diameter [mm]	front	200	200
	rear	600	600
use and storage temperature	–5°C to 40°C	–5°C to 40°C	–5°C to 40°C
static stability [°] – downhill, uphill, transverse	10°	10°	10°

ALUMINUM WHEELCHAIR – CRUISER LIBER

RF-20

Available colors:

○ white



- lightweight wheelchair from 10.6 kg
- swing away backrest
- adjustable gravity center
- height-adjustable footrest platform
- quick-release rear wheels

Intended use:

The Cruiser Liber wheelchair is intended for people with permanent disabilities limiting independent walking and maintained functionality in the upper limbs, enabling independent movement in a wheelchair.

TECHNICAL PARAMETERS							
reference number / name		CRUISER LIBER WHEELCHAIR – pneumatic wheels			CRUISER LIBER WHEELCHAIR – solid wheels		
seat width [mm]		420	450	480	420	450	480
total width [mm]		595	645	665	595	645	665
total length [mm]		820	820	820	820	820	820
total height [mm]		810	810	810	810	810	810
seat depth [mm]		370	370	370	370	370	370
backrest height [mm]		330	330	330	330	330	330
armrest height from the seat [mm]		220	220	220	220	220	220
max. safe load [kg]		100	100	100	100	100	100
basic wheelchair weight [kg]		11	11.4	11.5	10.6	11	11.2
static stability [°]	downhill	≥10°	≥10°	≥10°	≥10°	≥10°	≥10°
	uphill	≥10°	≥10°	≥10°	≥10°	≥10°	≥10°
	transverse	≥10°	≥10°	≥10°	≥10°	≥10°	≥10°
wheel diameter [mm]	front	102	102	102	102	102	102
	rear	595	595	595	595	595	595
ability to climb an obstacle [mm]		40	40	40	40	40	40
use and storage temperature		-10–45°C	-10–45°C	-10–45°C	-10–45°C	-10–45°C	-10–45°C

ALUMINUM RECLINING WHEELCHAIR – CRUISER UTILIS RF-11

Available
colors:

● gray



- swing away backrest
- swing away, height-adjustable and lift-up footrest
- detachable side panel
- anti-tip wheels
- flip-up, height-adjustable armrests
- comfortable headrest
- removable upholstery

Intended use:

The Cruiser Utilis wheelchair provides adequate support for the head and can therefore be used as a special wheelchair for the following conditions: paralysis and paresis, lower limb amputations, injuries, balance problems, muscular diseases, or other diseases limiting mobility.

It meets the needs of a user with locomotor system dysfunctions that prevent independent walking and requires stabilization of the body at least in a semi-recumbent and sitting position, as well as requiring a support system resulting from individual dysfunctions, given inability to maintain a sitting position on their own.

TECHNICAL PARAMETERS				
seat width [mm]		420	450	480
total width [mm]		630	660	690
total length [mm]		117	117	117
total height [mm]		129	129	129
seat depth [mm]		430	430	430
backrest height [mm]		490	490	490
armrest height from the seat [mm]		230	230	230
max. safe load [kg]		120	120	120
basic wheelchair weight [kg]		21.5	22	22.5
ability to climb an obstacle [mm]		60	60	60
wheel diameter [mm]	front	200	200	200
	rear	610	610	610
use and storage temperature		-10°~40°	-10°~40°	-10°~40°
static stability [°] – downhill, uphill, transverse		downhill ≥ 10°	downhill ≥ 10°	downhill ≥ 10°
		uphill ≥ 10°	uphill ≥ 10°	uphill ≥ 10°
		transverse ≥ 15°	transverse ≥ 15°	transverse ≥ 15°

STANDARD STEEL RECLINING WHEELCHAIR

YJ-011JA

Available
colors:

● gray



- swing away backrest
- swing away, height-adjustable and lift-up footrest
- detachable side panel
- headrest
- recessed gravity center, equipped with anti-tip elements
- foam cushion included
- soft height-adjustable headrest

Intended use:

A standard reclining wheelchair provides adequate head support and thus can be used as a special wheelchair for the following diseases: paralysis and paresis, lower limb amputations, injuries, balance problems, muscle diseases or other diseases limiting mobility.

It meets the needs of a user with locomotor system dysfunctions that prevent independent walking and requires stabilization of the body at least in a semi-recumbent and sitting position, as well as requiring a support system resulting from individual dysfunctions, given the inability to maintain a sitting position on their own.

TECHNICAL PARAMETERS

seat width [mm]		420	450	480
total width [mm]		620	645	670
total length [mm]		1510	1510	1510
total height [mm]		980	980	980
seat depth [mm]		450	450	450
backrest height [mm]		520	520	520
armrest height from the seat [mm]		430	465	485
max. safe load [kg]		120	120	120
basic wheelchair weight [kg]		25.30	25.40	25.80
ability to climb an obstacle [mm]		60	60	60
wheel diameter [mm]	front	200	200	200
	rear	600	600	600
use and storage temperature		-5°C–40°C	-5°C–40°C	-5°C–40°C
static axis [°] – downhill, uphill, transverse		10°	10°	10°

STANDARD ELECTRIC WHEELCHAIR

DY01109

Available
colors:

● blue



- magnesium alloy rear wheel hub
- automatic switching between manual and electric modes with a push button
- rear wheels adapted to ride on various surfaces: speed bumps, dirt road, stone surface, grass
- lithium battery
- 10 Ah battery
- range of 20 km
- multifunctional joystick with LCD display
- 360-degree easy-grip joystick
- speed adjustment
- real-time speed information
- Bluetooth and radio
- easy to fold
- compact dimensions after folding for easy transport
- wide adjustable hip belt
- foldable backrest
- lift-up footrests
- adjustable anti-tip wheels

Intended use:

This product is intended for the elderly, disabled and other groups of people with limited walking ability (except groups that cannot judge the environment they are in), suitable for indoor, courtyard, sidewalk and other places for short distances, it is an ideal means of transport for the elderly and the disabled.

TECHNICAL PARAMETERS

total length	99 cm
total width	68 cm
total height	91 cm
seat width	46 cm
front wheel diameter	8"
rear wheel diameter	20"
seat height	48 cm
seat depth	40 cm
backrest height	41 cm
width of the folded product	33 cm
product weight	25.5 kg
user weight	100 kg
max. speed	≤6 km/h
braking performance	≤1.5m
static stability	≥9°
dynamic stability	≥6°
min. turning angle	≤ 1200 mm
theoretical range	20 km/one battery

MOBILE LIFT – NEWTON

RF-1000

Available colors:

○ white




- adjustable width of the base spacing
- possibility to raise and lower the patient using a remote control
- emergency battery cut-off button
- rear wheel locks
- set with a standard M size seat included
- rotating hanger

Intended use:

The product provides support when moving patients from different places and medical devices. The lift is used to transport and move patients indoors. It helps while transporting people with total or partial disabilities, such as bedridden and elderly patients with reduced mobility. The correct use of the lift enables the transfer of elderly, disabled and paralyzed people, e.g. to a bed, toilet, bathroom, etc.

TECHNICAL PARAMETERS

dimensions	1200 × 595 × 1210 mm
case width	810 mm
battery (lead-acid) rated capacity	4.5 Ah
powered actuator	rated voltage 24 V
front wheel diam.	3 inches
rear wheel diam.	4 inches
user weight	150 kg
weight	42 kg
lifting height	1210-1820 mm
nominal voltage	24 V
pressure	8000 N
number of lifts with fully charged battery	80-120 times

CHOOSE ANY SEAT**Standard
RF-1010****Seat size**

A-46 cm | B-98 cm | C-95 cm

**With a head support
RF-1011****Seat size**

A-47 cm | B-88 cm | C-85 cm

**For bathing
RF-1012****Seat size**

A-46 cm | B-98 cm | C-95 cm

NOTES



ul. Staniewicka 14
03-310 Warszawa
Poland

info@rehafund.pl
www.rehafund.pl